

# ELISA PLATE ANALYSER

## MICRO READ 1000



CE

## **TABLE OF CONTENTS**

---

<b>1. GENERAL INFORMATION.....</b>	<b>5</b>
1.1. WARRANTY INFORMATION.....	5
1.2. TECHNICAL SERVICE .....	5
1.3. CONTACTS.....	5
<b>2. GENERAL SAFETY WARNINGS.....</b>	<b>6</b>
2.1. DANGER – WARNINGS SYMBOLS:.....	6
2.2. USE OF THE INSTRUMENT .....	7
<b>3. INTRODUCTION.....</b>	<b>8</b>
3.1. DESCRIPTION.....	8
3.2. SPECIAL FEATURES.....	8
3.3. TECHNICAL SPECIFICATION OF Micro read 1000: .....	10
<b>4. PACKAGING, TRANSPORT AND STORAGE.....</b>	<b>12</b>
4.1. GENERAL WARNINGS .....	12
4.2. PACKAGING .....	12
4.3. INSTRUMENT TRANSPORTATION.....	12
4.4. STORAGE OF INSTRUMENT.....	12
<b>5. INSTRUMENT DESCRIPTION.....</b>	<b>13</b>
5.1. INSTRUMENT WORKING PRINCIPLE .....	13
5.2. PERSPECTIVE VIEW .....	14
5.3. KEYPAD.....	16
<b>6. INSTALLATION PROCEDURE &amp; VERIFICATION CRITERIA.....</b>	<b>17</b>
6.1. UNPACKING INSTRUCTIONS .....	17
6.2. PLACING THE INSTRUMENT.....	17
6.3. POWER SUPPLY REQUIREMENTS.....	17
6.4. PROTECTIVE GROUNDING.....	17
6.5. PRINTER .....	17
6.6. START UP INSTRUCTIONS.....	19
6.7. TOUCH PANEL CHECK .....	20
6.8. MICRO TITER PLATE CARRIAGE .....	21
6.9. PLATE LOADING & PIPETTING PROCEDURE.....	21

---

6.10. READINGS CHECK .....	22
<b>7. PRECAUTIONS.....</b>	<b>23</b>
<b>8. PROGRAMMING MODES.....</b>	<b>24</b>
8.1. ABSORBANCE MODE .....	24
8.2. SINGLE STANDARD .....	27
8.3. CUT OFF MODE .....	28
8.4. REVERSE CUT OFF MODE.....	43
8.5. MULTI STANDARD MODE: .....	44
8.6. PERCENTAGE ABSORBANCE .....	55
8.7. UPTAKE .....	55
8.8. RATIO .....	56
<b>9. RECALLING AND RUNNING OF STORED TEST/PROGRAMS.....</b>	<b>57</b>
9.1. ABSORBANCE MODE: .....	58
9.2. CUT OFF MODE.....	65
9.2. A. INVALID ASSAY IN CUT OFF MODE .....	75
9.3. MULTI STANDARD MODE .....	78
9.3. A. INVALID ASSAY IN MULTI STANDARD MODE .....	92
<b>10. RERUNNING ACCEPTED TESTS/PROGRAMS.....</b>	<b>94</b>
<b>11. SAMPLE AND SAMPLE DUPLICATE.....</b>	<b>100</b>
<b>12. PATIENT IDENTITY / SAMPLE IDENTITY ENTRY.....</b>	<b>104</b>
<b>13. SAVE PLATE AND SAVE RESULTS.....</b>	<b>110</b>
<b>14. VIEW STORED DATA.....</b>	<b>114</b>
14.1. PLATE RESULT: .....	114
14.2. PATIENT RESULTS .....	118
<b>15. UTILITIES.....</b>	<b>120</b>
15.1. DATE & TIME SETTING.....	120
15.2. PROGRAM FILTER.....	121
15.3. PRINTER ON / OFF SETTING: .....	121
15.4. ENTER CLINIC NAME:.....	122
15.5. SHAKER (PLATE SHAKING MODE).....	122

---

15.6. COMMUNICATION SETTINGS:.....	125
15.7. CV DIAGNOSTICS: OFF .....	125
15.8. UPDATE FIRMAWARE: .....	125
15.9. BI-DIRECT .....	125
15.10. PASSWORD .....	126
<b>16. COMMUNICATION.....</b>	<b>129</b>
16.1. SEND RESULTS .....	130
16.2. SEND TEST RECORDS.....	132
16.3. RECEIVE TEST RECORDS .....	133
16.4. CHECK SERIAL COMMUNICATION.....	134
<b>17. TECHNICAL DIAGNOSIS.....</b>	<b>135</b>
17.1. CLEAR MEMORY .....	135
17.2. LAMP AMPLITUDE .....	136
17.3. PRINTER SELF TEST .....	136
<b>18. TROUBLESHOOTING.....</b>	<b>140</b>
<b>19. DECONTAMINATION.....</b>	<b>142</b>
19.1. DECONTAMINATION PROCEDURE .....	142
19.2. PURPOSE OF DECONTAMINATION .....	142
19.3. GENERAL CONSIDERATIONS .....	142
19.4. PROCEDURE.....	142
<b>20. SAFETY CLEARANCE CERTIFICATE.....</b>	<b>143</b>

## 1. GENERAL INFORMATION

### 1.1. Warranty Information

Each Instrument is completely tested and guaranteed for twelve months from delivery. The warranty applies to all the mechanical and electrical parts. It is valid only for proper installation, use and maintenance in compliance with the instructions given in this manual.

**Global Diagnostics B** will, at its discretion repair or replace parts, which may be found defective in the warranty period. The warranty does not include any responsibility for direct or indirect personal and material damages, caused by improper use or maintenance of the instrument.

Parts that are inherently subject to deterioration are excluded from the warranty. In case of defects due to misuse of the instrument, any travel and man-hour expenses will be charged extra.

In case of tenders warranty would be as per tender terms and conditions.

### 1.2. Technical Service

**Global Diagnostics B** is always accessible to the customers for any kind of information about installation, use, maintenance and others. When asking for service, please refer to this manual and report the data reported on the identification label (serial number).

Only qualified technicians are entitled to repair instruments. The user should carry out ordinary maintenance.

The technical service of **Global Diagnostics B** or an authorized service center with specialized technicians, with suitable instrumentation and original spare parts only, is always available for extraordinary maintenance (repair), under a yearly maintenance contract or on specific demand.

### 1.3. Disposal Instruction:

In case of removal or disposal of instrument, following instructions need to be followed

- Do not dispose in municipal waste; follow local regulations for instrument disposal.
- Plastic parts, Electronic PCBs and components can be recycled, so return back the instrument to manufacturer.

### 1.4. Contacts:

#### Manufacturer:

**Global Diagnostics B,**  
**Sijsjesstraat 4, 2440 Geel,**  
**Belgium.**  
**Tel: +32468220039**  
**info@globaldiagnosticsb.com**



## 2. GENERAL SAFETY WARNINGS

### 2.1. Danger – warnings symbols:

The following symbols are used to inform the user of the safety rules.



This symbol indicates generic danger. It means that, serious damage can occur to the operator if described precautions are not observed.



This symbol indicates HIGH ELECTRIC VOLTAGE. It is dangerous to touch any part having this label. Only qualified operators can access these components, after unplugging the instrument from the Supply.



This symbol indicates that the instrument involves the handling of samples, which can be infected (urine or human serum). In this condition, infection or contamination might occur. Pay attention to the general safety warnings when in presence of such biological substances. Use Protective clothes, gloves and glasses.



This symbol in the user manual indicates that damages to the instrument or erroneous results could occur if the given warnings are not followed.



This symbol indicates a portion, which is particularly important, and should be studied carefully.



This symbol indicates a Protective Earth or Ground terminal.

### General Symbols



Symbol for “Manufacturer”



Symbol for “ IN VITRO DIAGNOSTIC MEDICAL DEVICE”

## 2.2. Use of the Instrument

1. The instrument has to be used for the designed purposes under specified conditions, following proper procedures and safety rules by qualified personnel.
2. This manual contains instructions for operation performed by qualified personnel.
3. A qualified user has to make sure that environmental condition is suitable, the installation is correct, the use and maintenance are proper according to the general safety rules as well as to the particular precautions described in the manual. (However, he is not entitled to repair the instrument)
4. A qualified technician is entitled to maintain and repair the instrument using the original spare parts according to the given instructions.
5. Maintain room temperature and humidity as specified in the manual.
6. If the instrument is not used as described in the manual, the protection provided by the instrument may be impaired.
7. Alterations to the instrument are prohibited. The user is liable for any improper modification to the instrument, and for the deriving consequences.
8. Contact the MANUFACTURER service or authorized service center in case the instrument need extraordinary maintenance. Specialized technicians who will be able to repair the instrument using original spare parts will carry out the maintenance.
9. This IVD equipment complies with the emission and immunity requirements as per IEC61326 series.



10. **Warning:** This equipment has been designed and tested to CISPER11 Class A. In a domestic environment it may cause radio interference, in which case, you may need to take measures to mitigate the interference."

11. An advisory that the electromagnetic environment should be evaluated prior to operation of the device.



**Warning:** Do not use this device in close proximity to sources of strong electromagnetic radiation (e.g. unshielded international RF sources), as these may interfere with the proper operation.

---

## 3. INTRODUCTION

### 3.1. Description

**MICRO READ 1000** is a user friendly micro plate Analyser. It is compact & lightweight. It is designed to measure and interpret enzyme immunoassay results, both monochromatically and bichromatically. It is intended for in vitro diagnostic use.

### 3.2. Special Features

The **MICRO READ 1000** can accommodate a flat bottom as well as a round configuration. The carriage is designed in a way that the plate automatically moves smoothly and positions itself accurately in the optical measurement path. Readings are taken continuously. The average value is calculated and results are presented according to the option selected.

The **MICRO READ 1000** operates on a WIDE voltage (90-270 volt). This eliminates the need for an external voltage stabilizer.

The **MICRO READ 1000** has a special provision, which allows it to be used even when a printer is not available. Readings can be conveniently noted down manually.

The other special features of **MICRO READ 1000** are as follows:

- Option of Lamp saving mode.
- Selection of both primary and secondary filters.
- Latest technology with battery backup for 250 tests, and 100 plate data.
- Robust built in 52-column thermal printer with 384 stationary heads.
- Unique circuitry for long lamp life.
- Alphanumeric Patients ID entry.
- Editing of saved tests.
- Human machine user interface: Touch panel, Keypad
- Multi-standard curve up to 12 standard calibrations with one blank optional.
- Access to test by touch of key.
- It takes 16 secs to read the plate.
- With help of color combination, we can distinguish the calibrator, control, high and low sample
- Blank is optional.
- Setting of the Date and Time.
- Capable of storing, deleting and recalling tests.
- Multiple calibrator modes.
- Selection of duplicates for both calibrators and samples.



- Extensive software for Cut Off mode.
- Selection of Positive, Equivocal, Negative cut-off.
- Several pre-programmed calculation modes will help to facilitate data processing of enzyme immunoassays. These are menu driven modes for simple and error free operations.

- ✓ ABSORBANCE MODE
- ✓ SINGLE STANDARD MODE
- ✓ CUT-OFF MODE
- ✓ MULTISTANDARD MODE
- ✓ % ABSORBANCE MODE
- ✓ UPTAKE
- ✓ RATIO

**3.3. Technical Specification of MICRO READ 1000:**

Human Machine Interface	TOUCH PANEL / KEYPAD
Linear measurement range	0.001 to 3.000 Absorbance Units (A)
Photometric Accuracy	± 2% or 0.007 whichever is higher, from 0 to 1.5 A ± 3% from 1.5 A to 3.0 A
Drift	<0.005 A/hr
Photometric Linearity	2.5 A
Optical measurement	2 Channel
Filters	
Type of filter	Narrow band Interference
Wave Length	405nm, 450nm, 492nm,560 nm 630nm, with two optional filters namely – XXX and YYY (editable)
Half Bandwidth	10nm ± 2nm
Selection	Automatic by Stepper Motor
Light Source	Tungsten halogen lamp, 20 Watts
Display	7" TFT
Curve Plotting	STN 480*800 pixels Graphical Representation on Printer
Plate Carrier Movement	Precisely through the stepper motor
Printer	Built in Thermal Printer 52 columns
Memory	2MB FRAM
Analysis Mode	Absorbance Single Standard Cut-off Multi-Standard % Absorbance Uptake Ratio

Connectivity / RS232 Serial Port / USB	9600 baud, 8 data, 1 stop, no parity bits / USB
Power Wattage Voltage	75 Watts 115-230 Volt $\pm$ 10%, 50/60 Hz
Operation Position	On horizontal flat, rigid and vibration free surface
Operating Conditions Temperature Relative Humidity	From + 18°C to 35°C Up to 80%
Storage Conditions Temperature Relative Humidity	From -10°C to 40°C Up to 80%
Enclosure	ABS Fire Retardant
Size (cm)	36 x 36 x 22 (lxbxh)
Well Type	Capable of reading U,V and Flat type wells
Weight (Approx.)	10 Kg

## **4. PACKAGING, TRANSPORT AND STORAGE**

### **4.1. General Warnings**

Instrument has to be decontaminated before packing for transportation.

### **4.2. Packaging**

Packaging is needed whenever the instrument is to be transported or shipped by courier or other purposes.

 To pack the instrument the following instructions has to be followed:

1. Decontaminate the instrument as explained on decontamination chapter of this manual.
2. Put the instrument into the original packaging box; Instrument has to be properly protected by plastic protective material. Put copy of Safety clearance certificate (copy of Safety clearance certificate is attached at the end of this manual)
3. Mark the packaging with address, instrument identification and warning labels

### **4.3. Instrument Transportation**

The transportation of the instrument in unpacked condition must be limited within the room where it is used, to avoid damage.

### **4.4. Storage of Instrument**


Before storing the instrument for a long period, pack it carefully as described above and store indoors.

Relative humidity has to be less than 85%, and temperature between -10°C and 40°C.

## 5. INSTRUMENT DESCRIPTION

### 5.1. Instrument Working Principle:


- The following diagram represents the main functional elements of the instrument:
- White light produced by the lamps is focused into a beam by the lens & passes through the sample. Part of the light is absorbed by the sample & the remaining light is transmitted. It is filtered by interference filters & focused onto the photodiodes. The photodiode converts the received light in to an electrical signal which is in-turn transformed into digital form, from which the microprocessor calculates the absorbance, taking in account of the blank & Bichromatic selection.




---

## 5.2. Perspective View


### Front View:



**Rear View:**

- 1) Serial RS232 output
- 2) Cooling Fan
- 3) USB output
- 4) ON/OFF switch
- 5) SMPS

### 5.3.Keypad



- 1) 'Yes' key is used to select any 'YES' option on screen directly.
- 2) 'No' key is used to select any 'NO' option on screen directly.
- 3) 'Print' key is used to get a printout of current screen displayed.
- 4) 'Feed' key is used to advance the printer paper by 1 line.
- 5) Navigation keys are used to select any option available on current screen.
- 6) 'ENTER' key is used to run any option selected by navigation keys.
- 7) 'ESC' key is used for escaping from any screen.



## 6. INSTALLATION PROCEDURE & VERIFICATION CRITERIA

### 6.1. Unpacking Instructions

Check accessories as per packing list.

Kindly store all packaging materials so as to use it to repack and ship for maintenance or servicing.

### 6.2. Placing the Instrument

The instrument has to be placed on a level bench.

Room temperature has to be between 10 to 35°C with a relative humidity below 85%.

Protect the instrument from direct sunshine.

### 6.3. Power Supply Requirements

Plug the instrument into a power source by the locally available approved plug in cable.

Power cord should be CE, CSA, and UL marked.

**115 - 230 Volt  $\pm$  10%, 50-60 Hz**

### 6.4. Protective Grounding




Please make sure that electrical power source is properly grounded.

### 6.5. Printer

“micro read 1000”, an ELISA Analyser is equipped with a built in 52 columns easy load Thermal printer. Procedure to load the paper is as follows:

#### Opening the Paper Cover:

Slide the LEVER towards back of the printer to open lock of the Paper Cover.




**Closing the Paper Cover + replacement of paper roller:**

Open the Paper Flap .

Remove the paper roller from the slot.

Now place the Thermal side of the paper roll at the top into the slot provided and then close the paper flap in the direction as shown in the following figure, till you get the locking sound. Use Thumb impressions to push the cover .

Press the paper feed switch until the paper feeds straight and smoothly.

**HINT:**

When the paper is set correctly and when the closing of the paper cover is proper the FEED LED will not glow.


It keeps on flashing when the data is being printed.

When the paper is almost finished, red lines appear on both sides of the paper.


## 6.6. Start up Instructions

Switch on the instrument.

The instrument initializes all the parameters internally and carries out a power on self-test. During Initializing it displays the ADDRESS screen. After address screen it will display the string: "Please wait... system initialization".



The instrument will print the model name "micro read 1000", Version number, Clinic Name, current Date with Time.



User can select a saved test by first touching the required "Test Name" on the screen and then touch the "VIEW" option to carry out further operation on the test.

"RUN" option is used to directly run the selected test.

"DEL": Deletes the selected test.

If the numbers of tests are more than 18, "NXT" and "PRV" options can be used for browsing through the list pages.

"PRN": For printing the test screen.

"MENU" used to select Main menu screen.

"NEW" this option helps the user to create new test

## 6.7. Touch Panel Check

Micro read 1000 provides a *TFT LCD panel* and a *KEYPAD* for easy user interface.

The Menus are displayed and the text box of the parameter forms the TOUCH ZONE.

### Touch screen Layout

NAME	MODE COFF	PRI 405	SEC 0	ESC
BLK N 0 0.000	NC N 0 0.000	PC N 0 0.000	LC N 0 0.000	SAVE
CC N 0 0.000	QC 0			
CUT-OFF INDEX 0.000 0.000		GREY ZONE 0.000 0.000		
CUT-OFF EQUATION				

Above is the generic representation of a Test Screen. The Highlighted zones are TOUCH ZONES, which are active. On touching the “Touch Zone” of a parameter, a sub menu/menu is displayed or the requested action is carried out. The rest of the “Touch Zone” is deactivated.

**For Example:** To activate the selection


Enter the primary filter value to touch any point in the shaded area “PRI” on the screen.

On proper selection the analyser responds with blinking of the parameter text and also the TOUCH ZONE and a submenu is displayed.

To enter Test Name: Touching the “NAME” touch zone provides an alphanumeric screen.

Enter the Test name by touching the Touch zone of that variable. The selected value blinks and is displayed next to the parameter.

## 6.8. Micro Titer Plate Carriage



The instrument is provided with a micro titer plate carriage to move the micro titer plate inside. This carriage is driven by stepper motor with timing belt. It places each well of plate exactly below the optical path of each channel. The plate carrier is moved by a well controlled stepper motor drive.

## 6.9. Plate Loading & Pipetting Procedure


### Controls and Sample pipetting procedure in case of CUT OFF:

Maximum number of Blanks = 5

Maximum number of Controls:

1. Negative Control = 05
2. Positive Control = 05
3. Low Positive Control = 05

Pipetting sequence should be as shown below.



## 6.10. Readings Check

Checking of readings should be done through controls reading within range specified in data sheet of controls (care should be taken while preparing and pipetting controls and reagents, reagent and control expiry dates need to be checked).

## 7. PRECAUTIONS



1. Do not use any sharp objects on the Touch Screen. Always use the STYLUS provided to operate the touch panel.



2. Always check for proper grounding before installation. Never operate the instrument when ground wire is removed.

3. Do not attempt to open the instrument and make repairs without proper technical training. Do not allow unauthorized persons to operate or repair the instrument.

4. Use a clean plate and follow the instructions for blanking and standardizing. Do not read any wells containing bubbles or dust particles.

5. The volume of sample, calibrators and blanks should be identical for correct readings. The absorbance is proportionate to the path length. Pipetting should also be proper.

6. Monitoring of the printed values or displayed values during operation may help detect an error. Check the linearity and calibration of the instrument regularly against some standard reference.

7. Check the micro wells before use. They should be scratch-free. The micro well track in the **micro read 1000** has been designed in such a way that the micro wells are totally protected from scratches.



8. Recheck the reading of high OD (above 2A).

9. Place Plate carefully on the tray.

10. Ensure that the main power switch is in OFF position before connecting.

11. Plug the instrument to the AC mains. Confirm proper grounding for trouble free operation.

12. Connect the printer only when the instrument is OFF.


## 8. PROGRAMMING MODES

### 8.1. Absorbance Mode


In this mode the instrument gives only absorbance values of all wells of the respective plate.

#### Programming a New Test:


##### A. Select "NEW"



##### B. Select "MODE"



##### C. Select mode of operation "Absorbance"





**D. Enter “Test Name”**


The diagram illustrates the process of entering a test name. On the left, a menu screen displays 'NAME', 'MODE ABS', 'PRI 405', 'SEC 0', 'ESC', 'BLK N 0.000', and 'SAVE'. An arrow points to the right, where a keyboard interface is shown. The keyboard has a grid of letters (A-Z, numbers 0-9, and symbols like ., -, +, \_) and function keys (ENT, CLR, SPC). The 'NAME:' field is populated with 'ABSORB'.

**E. Select primary and secondary filter**

The diagram illustrates the process of selecting a filter. On the left, the same menu screen as in section D is shown. An arrow points to the right, where a filter selection screen is displayed. The screen shows 'NAME', 'MODE ABS', 'PRI 405', 'SEC 0', 'ESC', 'BLK N 0.000', and 'SAVE'. Below these fields, the number '405' is displayed, and a grid of filter options is shown: 000, 405, 450, 492, 630, XXX, and YYY.


Select respective filter value which is display in above screen.

**F. If blank is required select “BLK”**



If you select “BLANK Y”, it will read first strip first well ‘A1’ as a blank and subtract the absorbance of blank well from all other well’s absorbance.

After entering all parameter click on save button.



## 8.2. Single Standard

In this mode the instrument accepts the calibrator singly or in duplicate and then calculates the concentration based on the single point standard curve passing through the point 0.000

A single calibrator/standard of a known concentration is used to calibrate the instrument so that the concentration of unknown samples can be calculated according to Beer's Law. The absorbance are read at user selected wavelengths. If Blank is selected, instrument will automatically blank on the first well and subtract its absorbance from each subsequent well. The second well is treated as the calibrator/standard well. The third well is also treated as calibrator/standard well if the calibrator/standard is in duplicate.

### CALCULATION:

Sample Concentration = (Calibrator concentration /Calibrator Abs.) x Sample Abs

NAME SIN_STD	MODE SSTD	PRI 450	SEC 0	ESC
BLK Y 0.000	CAL 1.000	HI CTRL 0	LO CTRL 0	SAVE
NORMAL RANGE				

Entry of all the parameters is similar to Multi standard mode (Please refer Multi standard for entry of parameters except unit selection )

### 8.3. Cut Off Mode

In this mode Cut-Off point is determined for interpretation of specimens as per formula given in the reagent manual. The negative controls are read followed by the positive controls, Cut-off control, low positive control. Blanking on the first well is optional. The instrument calculates the average of the negative controls and the average of the positive control, Cut-off controls and Low positives are also calculated.

#### Programming a New Test:

##### A. Select "NEW"

SIN
TSH
UPT
CPC
KIN
CRP
NXT
PRV
NEW
VIEW
DEL
MENU
RUN
PRN

##### B. Select "MODE"

NAME	MODE	PRI	SEC	ESC
	ABS	405	0	
BLK				SAVE
N				
0.000				

##### C. Select mode of operation "CUT-OFF"

NAME	MODE	PRI	SEC	ESC
ABSORB	ABS	405	0	
BLK				SAVE
N				
0.000				
MODE				
ABSORBANCE	SINGLE-STANDARD			
ABS	SSTD			
CUT-OFF	MULTI-STANDARD			
COFF	MSTD			
PERCENT ABS	CUT-OFF REVERSE			
%ABS	CUTREV			
UP-TAKE	RATIO			
UPT				

**D. Entry of test name**

Select touch Zone “Name:”


NAME	MODE COFF	PRI 405	SEC 0	ESC
BLK N 0 0.000	NC N 0 0.000	PC N 0 0.000	LC N 0 0.000	SAVE
CC N 0 0.000	QC 0			
CUT-OFF INDEX 0.000 0.000		GREY ZONE 0.000 0.000		
CUT-OFF EQUATION				

Alphanumerical characters will appear on the screen as shown below after touching “Name” on the screen. Enter the name of the test by selecting the individual characters one by one. After selecting all the characters select Enter option in touch screen for confirmation of test Name.

NAME	MODE COFF	PRI 405	SEC 0	ESC				
BLK N 0 0.000	NC N 0 0.000	PC N 0 0.000	LC N 0 0.000	SAVE				
CC	QC							
NAME: TSH								
A	B	C	D	E	F	G	H	ENT
I	J	K	L	M	N	O	P	CLR
Q	R	S	T	U	V	W	X	SPC
Y	Z	.	_	-	+	1	2	
3	4	5	6	7	8	9	0	

**E. Selection of Primary and secondary filters.**

Select zone “PRI” on the touch screen




Select the filter value for primary filter as per details given in the reagent manual and after selection of “PRI” filter it will come back to the test parameters screen.

**NOTE:** Same way you can select the filter Values for secondary filter.


Select the required filters either monochromatic or bichromatic by selecting the filter values for primary and secondary filter.

**F. Entry of Blank and Blank QC values**

Select “BLK” to enter blank details.




**G. After selecting “BLK”, enter total number of blanks by selecting “NO OF BLANK”.**  
**(Maximum 5 blank).**




**H. Enter the QC value of the blank by selecting the touch Zone “QC”**

For entering QC value select “QC N”. After selecting QC select “VALUE OF BLANK” to enter the QC value of blank.




**I. Entry of Negative control and its QC values**


Select touch Zone “NC”




**J. After selecting “NC”, below screen will appear on the display.**



**Enter number of Negative Controls by Selecting Zone “NO OF NC” (Maximum 5).**



**Enter the QC value of the Negative Control by selecting the touch Zone “QC N” and enter the value of Y given in the reagent manual by selecting Zone “VALUES OF NC”.**





**K. After the entry of “No. of NC’s”, select ESC. The below screen will appear on the display.**


NAME TSH	MODE COFF	PRI 450	SEC 630	ESC
BLK Y 1 0.100	NC Y 3 0.150	PC N 0 0.000	LC N 0 0.000	SAVE
CC N 0 0.000	QC 0			
CUT-OFF INDEX 0.000 0.000		GREY ZONE 0.000 0.000		
CUT-OFF EQUATION				

**L. Entry of Positive Controls (PC).**


Select touch Zone “PC”.

NAME TSH	MODE COFF	PRI 450	SEC 630	ESC
BLK Y 1 0.100	NC Y 3 0.150	PC N 0 0.000	LC N 0 0.000	SAVE
CC N 0 0.000	QC 0			
CUT-OFF INDEX 0.000 0.000		GREY ZONE 0.000 0.000		
CUT-OFF EQUATION				


**M. After selecting “PC”, below screen will appear on the display.**



**Enter number of Positive Controls by Selecting Zone “NO OF PC” (Maximum 5).**



**Enter the QC value of the Positive Control by selecting the touch Zone “QC N” and enter the value of Y given in the reagent manual by selecting Zone “VALUES OF PC”.**



N. After entering the details of PCs, escape from the screen by selecting “ESC” option.

Below screen will appear on the display:

NAME TSH	MODE COFF	PRI 450	SEC 630	ESC
BLK Y 1 0.100	NC Y 3 0.150	PC N 0 0.000	LC N 0 0.000	SAVE
CC N 0 0.000	QC 0			
CUT-OFF INDEX 0.000 0.000		GREY ZONE 0.000 0.000		
CUT-OFF EQUATION				

**NOTE:** Same way you can enter values for LC (lower positive controls) and CC (Cut of Controls)

**O. Entry of Cut off absorbance equation:**

Select touch Zone “CUT-OFF EQUATION”

NAME TSH	MODE COFF	PRI 450	SEC 630	ESC
BLK Y 1 0.100	NC Y 3 0.150	PC N 2 0.500	LC N 0 0.000	SAVE
CC N 0 0.000	QC 0			
CUT-OFF INDEX 0.000 0.000		GREY ZONE 0.000 0.000		
CUT-OFF EQUATION				

After selecting the “CUTABS” below screen will come on the display.

NAME TSH	MODE COFF	PRI 450	SEC 630	ESC
BLK Y 1 0.100	NC Y 3 0.150	PC Y 2 0.500	LC N 0 0.000	SAVE
CC N 0 0.000	QC 0			
CUT-OFF INDEX 0.000 0.000		GREY ZONE 0.000 0.000		
CUT-OFF EQUATION				

Now enter the equation as given in the reagent manual or else select ENT option in order to come out of the CUTABS screen.

Whereas, NC = Mean of NCs

PC = Mean of PCs

LC = Mean of LCs

CC = Mean of CCs

ALg = Anti LOG

Sqrt = Square root

Lg = LOG

Ln = Natural LOG

Abs. = any absolute value


ALN = Anti Natural LOG


**P. QC1, QC2, QC3 & QC4. can be used for any QC Checking of the test provided in reagent manual.**

For example some kit may give QC check equation like  $|PCx - NCx| > 0.500$ ,  $PC / NC > 15$ , etc. This equation can be entered by selecting these QC options.

NAME TSH	MODE COFF	PRI 450	SEC 630	ESC
BLK Y 1 0.100	NC Y 3 0.150	PC Y 2 0.500	LC N 0 0.000	SAVE
CC N 0 0.000	QC 0			
CUT-OFF INDEX 0.000 0.000		GREY ZONE 0.000 0.000		
CUT-OFF EQUATION NC+PC				

Select touch zone 'N'. It will change to 'Y', as shown in the below screen. Now select touch zone 'QC1' and enter QC equation.





Same way you can enter QC equations for QC2, QC3 and QC4.

**\*NOTE:**

On screen it is mentioned that “QC1 > Y”. In most of the reagent inserts it is mentioned that QC value should be greater than some constant value. Sometimes in the box insert it is mentioned that the value should be less than some constant value. In such case, you can change the symbol from “QC1 > Y” to “QC1 < Y”, only by selecting a symbol of comparison.

If any QC check is selected to ‘YES’ option, the system will not save the test without entering the QC equation and QC check value of equation.

**Q. Interpretation of samples**

Enter RANGE for interpretation of results either in percentage (Greyzone) or positive and negative values of INDEX, as a ratio of (Sample Abs. / Cutoff Abs.). If Cutoff Index range or Greyzone is not mentioned do not enter any values. Instrument will do the interpretation with references to the Cut off absorbance.

**Q.1. Interpretation by GREYZONE (percentage/range).**

Select touch zone 'GREYZONE

Select touch zone 'PERCENTAGE'

/ 'RANGE'

NAME TSH	MODE COFF	PRI 450	SEC 630	ESC
BLK Y 1 0.100	NC Y 3 0.150	PC Y 2 0.500	LC N 0 0.000	SAVE
CC N 0 0.000	QC 1			
CUT-OFF INDEX 0.000 0.000		GREY ZONE 0.000 0.000		
CUT-OFF EQUATION NC+PC				

GREY ZONE	ESC
PERCENTAGE 0.000	
RANGE 0.000-0.0000	

GREY ZONE	ESC
PERCENTAGE 0.000	
RANGE 0.000-0.0000	
ENTER NUM: 10	
1 2 3 4 5 6 ENT	
7 8 9 0 . CLR	

Enter Greyzone PERCENTAGE value  
Select ESC option after entering Greyzone % value to come out to the parameter screen.

RANGE	ESC
HIGH VALUE 0.000	
LOW VALUE 0.000	
ENTER NUM: 10	
1 2 3 4 5 6 ENT	
7 8 9 0 . CLR	

Enter Greyzone RANGE value  
Select ESC option after entering Greyzone range value to come out to the parameter screen.

- 
- NOTE: 1. If you don't select any option for interpretation of sample results, the instrument will take Cutoff Absorbance as a reference.**
- 2. when you are using a 'Greyzone' option the 'Cut off Index Range' option will not work.**


It will give "Positive" remark to sample absorbance greater than cutoff and "Negative" to sample absorbance less than cutoff.

NAME TSH	MODE COFF	PRI 450	SEC 630	ESC
BLK Y 1 0.100	NC Y 3 0.150	PC Y 2 0.500	LC N 0 0.000	SAVE
CC N 0 0.000	QC 1			
CUT-OFF INDEX 0.000 0.000		GREY ZONE 10.00		
CUT-OFF EQUATION NC+PC				



**Q.2. Interpretation by ‘Cutoff Index Range’**

It is also provided to do the interpretation using ‘Cutoff Index’ by entering normal range for ‘Positive’ and ‘Negative’.



In case of ‘Cutoff Index Range ’,  
 It gives the ‘Positive’ remarks to sample having index value, a ratio of ...  
 (Sample abs. / Cutoff Abs.) greater than or equal to entered ‘POSITIVE’ value and  
 it gives ‘Negative’ remarks to sample having index value, a ratio of  
 (Sample abs. / Cutoff Abs.) less than entered ‘NEGATIVE’ value.  
 The sample having index value in between ‘POSITIVE’ and ‘NEGATIVE’ range get remarks ‘Equivocal’ =  
 ‘EQU’.

After entering the Interpretation values (ie. conditions for interpretation) the screen will be displayed as follows -

NAME TSH	MODE COFF	PRI 450	SEC 630	ESC
BLK Y 1 0.100	NC Y 3 0.150	PC Y 2 0.500	LC N 0 0.000	SAVE
CC N 0 0.000	QC 1			
CUT-OFF INDEX 0.000 0.900		GREY ZONE 0.000 0.000		
CUT-OFF EQUATION NC+PC				

For example:

If you enter POS  $\geq 1.0$  and NEG  $\leq 0.9$ .

Then, the sample will get POSITIVE remarks having INDEX Value, (Sam. Abs./ Cutoff Abs.)  $\geq 1.000$ .

The sample will get NEGATIVE remarks having INDEX Value, (Sam. Abs./ Cutoff Abs.)  $\leq 0.900$ .

The sample will get EQUIVOCAL remarks having INDEX Value, (Sam. Abs./ Cutoff Abs.) in the range from 0.900 to 1.000.

**R. To save the test after entering all the parameters, select ‘SAVE’ option.**

NAME TSH	MODE COFF	PRI 450	SEC 630	ESC
BLK Y 1 0.100	NC Y 3 0.150	PC Y 2 0.500	LC N 0 0.000	SAVE
CC N 0 0.000	QC 1			
CUT-OFF INDEX 0.000 0.900		GREY ZONE 0.000 0.000		
CUT-OFF EQUATION NC+PC				

## 8.4. Reverse Cut Off Mode

For Reverse Cut Off mode, select touch zone 'CUT-OFF REVERSE' key.

NAME	MODE CUTREV	PRI 405	SEC 0	ESC
BLK N 0 0.000	NC N 0 0.000	PC N 0 0.000	LC N 0 0.000	SAVE
CC N 0 0.000	QC 0			
CUT-OFF INDEX 0.000 0.000		GREY ZONE 0.000 0.000		
CUT-OFF EQUATION				

In Reverse Cut Off Mode, programming of test is the same as Positive Cut off Mode. Only QC check conditions for blanks and controls get reversed. In the 'INTERPRETATION', sample absorbance, which is lesser than cut off absorbance, gets remarked POSITIVE. Sample absorbance higher than cut off absorbance, gets remarked NEGATIVE. This is exactly the opposite to normal Cut Off Mode. Same way in case you select 'Cutoff Index Range'.

**.5. Multi Standard Mode:**

In this mode the instrument accepts a maximum of 12 calibrators and calculates concentration based on the best-fit curve. Graph is printed with Absorbance on Y-axis and concentration on X-axis.

**Programming a New Test:**

**A. Select “NEW”**

**B. Select “MODE”**

**C. Select mode of operation “MULTI-STANDARD”**

---

After selecting Multi standard mode it will display following screen –

The screenshot displays a control panel with the following elements:

NAME	MODE MSTD	PRI 405	SEC 0	ESC
BLK N 0.000	CAL 2	HI CTRL 0	LO CTRL 0	SAVE
NORMAL RANGE		UNIT		

#### D. Entry of Test Name

Test Name can be entered by Selecting Touch Zone “Name:”

Procedure for entering the test name is same as that of the Cut Off mode. Refer Test Name entry in Cut off mode for further details.


#### E. Selection of Primary and secondary filters

Primary and secondary filters can be selected by Selecting Touch Zone of “PRI” and “SEC”


Procedure for Selecting filter is same as that of the Cut Off mode. Refer selection of primary and secondary filters in Cut off mode for further details.

### F. Selection of Blank and duplicate blank

Select Touch Zone “BLK”. It will display below screen.




After selection of blank and duplicate blank it will display following screen.



### G. Entry of Standard Concentration

Select touch Zone “CAL” to enter Number of calibrators and its concentration.




After selecting “CAL”, it will display following screen:

Select Touch Zone “NO OF CAL”. Numeric screen will be displayed at the bottom of the screen.

Select the number of calibrators (User can select min 2 and max. 12 standard in this mode).


Whereas in case of 4 Parameter, user can select min. 4 and max.12 calibrators.



### Concentration Entry:

After entering the number of calibrators, concentration and absorbance column will come on the display as shown below.

To enter concentration values select touch zone of CON 1 in Conc. column corresponding to the standard and enter the concentration value. (Note: 1. The concentration values should be either in ascending or descending order).



### Graph type along with its scales:

There are five types of graphs in Multistandard. They are mentioned as follows -

- 1 . Linear (LINEAR)
- 2 . Point to Point (PT TO PT)
- 3 . 4 Parameter (4 PARAM)
- 4 . polynom
- 5 . CSPLINE

There are five types of scales for X and Y axis

LOGABS vs CONC (X-axis = Concentration v/s Y-axis = LOG of ABS.)

ABS vs LOGCONC (X-axis = LOG of concentration v/s Y-axis = Absorbance)

LOGITABS VS LOGCONC (X-axis = LOG of concentration v/s Y-axis = LOGIT of Absorbance)

ABS vs CONC (X-axis = Concentration v/s Y-axis = Absorbance)

LOGABS vs LOGCONC (X-axis = LOG of concentration v/s Y-axis = LOG of absorbance)



**Selection of graph:**

Select graph type as per the details given in the reagent manual by selecting Touch Zone “GRAPH”. After Selecting the “Graph” option present in the calibrator parameter screen, it will display following screen.

NO OF CAL 7				PLOT	ESC
CON 1 0.000	ABS 1 0.000	CON 7 0.000	ABS 7 0.000	SCALE ABS VS CON	
CON 2 0.000	ABS 2 0.000			GRAPH PT TO PT	DUP N
CON 3 0.000	ABS 3 0.000				QC N
CON 4 0.000	ABS 4 0.000			GRAPH ESC	
CON 5 0.000	ABS 5 0.000			PT TO PT	
CON 6 0.000	ABS 6 0.000			LINEAR	
				4 PARAM	
				CSPLINE	


**Selection of scale:**

Select scale type as per the details given in the reagent manual by selecting Touch Zone “SCALE”

NO OF CAL 7				PLOT	ESC
CON 1 0.000	ABS 1 0.000	CON 7 0.000	ABS 7 0.000	SCALE ABS VS CON	
CON 2 0.000	ABS 2 0.000			GRAPH PT TO PT	DUP N
CON 3 0.000	ABS 3 0.000				QC N
CON 4 0.000	ABS 4 0.000			SCALE ESC	
CON 5 0.000	ABS 5 0.000			ABS VS CON	
CON 6 0.000	ABS 6 0.000			ABS VS LOGCON	
				LOGABS VS CON	
				LOGABS VS LOGCON	
				LOGIT VS LOGCON	


**Selection of Duplicate Calibrators:**

To select Duplicate Calibrator, select Touch Zone “DUP” present in calibrators (“CAL”) parameter screen. After selecting “DUP” it will toggle to yes “Y” for Duplication of calibrators.




**Entry of QC:**


Select QC to enter the qc value.



Select touch zone ‘N’. It will display below screen. Now select touch zone ‘QC VAL’ and enter QC value.




Screen after entry of calibration and selection of duplicate calibration, graph and scale will be displayed as follows -




**H. Selection of Controls**

For High control select touch zone 'HI CTRL'. And enter the number of controls by selecting touch zone "NO OF CONTROLS"




Enter range of control “HIGH” and “LOW” as per reagent manual.



After entering HIGH and LOW value click on ESC to come back in test parameter screen.

Similarly, for Low Control select touch zone ‘LO CTRL’.



(Note: Entry of LO CTRL is similar to that of HI CTRL)

**I. Normal Range selection for interpretation of samples:**

For interpretation of specimen result (HI or LO), enter normal range given in REAGENT MANUAL. Select “NORMAL RANGE”. It will display “NORMAL RANGE” SCREEN. Then select “HIGH>” and “LOW<” to enter normal range.

NAME MUL_4P	MODE MSTD	PRI 450	SEC 0	ESC
BLK DUP Y 0.000	CAL DUP 7	HI CTRL 1 25.03 10.54	LO CTRL 1 16.61 5.45	SAVE
NORMALRANGE		UNIT		

All the results whose concentration is coming above the High Range will be considered as HI and coming below the Low Range will be consider as LO. The results having concentration in between High and Low range will get remarks 'Equivocal' = 'EQU'.

Enter Normal Range for positive and negative samples.

NORMAL RANGE	ESC
HIGH VALUE 0.000	
LOW VALUE 0.000	

→

NORMAL RANGE	ESC					
HIGH VALUE 15.25						
LOW VALUE 9.46						
ENTER NUM:						
1	2	3	4	5	6	ENT
7	8	9	0	.	CLR	

### J. Selection of unit

This option is used to select the unit for the result.

NAME MUL_4P	MODE MSTD	PRI 450	SEC 0	ESC
BLK DUP Y 0.000	CAL DUP 7	HI CTRL 1 25.03 10.54	LO CTRL 1 16.61 5.45	SAVE
NORMALRANGE 15.25 9.46		UNIT		

UNIT			ESC
No Unit	mg/ml	ng/ml	mIU/ml
pg/ml	ug/ml	UM/L	nm/L
U/ML	uIU/ML	ug/dL	g/dl
ng/dL	ppb	NEW	

### K. After entering all the parameters save the test by selecting “SAVE” option.

NAME MUL_4P	MODE MSTD	PRI 450	SEC 0	ESC
BLK DUP Y 0.000	CAL DUP 7	HI CTRL 1 25.03 10.54	LO CTRL 1 16.61 5.45	SAVE
NORMALRANGE 15.25 9.46		UNIT Mg/ml		

(\*Same way you can make a test for Single standard, % Absorbance and uptake. EXCEPT selection of UNIT)

## 8.6. Percentage Absorbance

The Percentage Absorbance Mode requires one calibrator (read singly or duplicate). In this mode, calibrator is considered to have a concentration of 100%. The absorbance's of unknown samples are read and compared to the calibrator absorbance, and reported as % concentration of calibrator.

Refer 8.5 for parameter entry.

NAME	MODE %ABS	PRI 405	SEC 0	ESC
BLK N 0.000	CAL 100.0	HI CTRL 0	LO CTRL 0	SAVE
NORMAL RANGE				

## 8.7. Uptake

In this mode the instrument accepts the calibrator singly or in duplicate and then calculates the concentration based on the single point standard curve passing through the point 0.000

A single calibrator/standard of a known concentration is used to calibrate the instrument so that the concentration of unknown samples can be calculated according to Beer's Law

Sample Concentration =  $\frac{\text{Calibrator Absorbance} * \text{Calibrator Concentration}}{\text{Sample Absorbance}}$

Sample Absorbance


NAME	MODE UPT	PRI 405	SEC 0	ESC
BLK N 0.000	CAL 0.000	HI CTRL 0 0.000	LO CTRL 0	SAVE
NORMAL RANGE				

Please refer 8.5 for parameter entry.

## 8.8. RATIO

The Ratio Mode is similar to multi standard mode and one more parameter is present in test parameter screen i.e. Ratio option present in test parameter screen. There are 2 types of Ratio

1. NORMAL: Instrument gives the absorbance of each and every well (A/A0).
2. %: Instrument gives the percentage of absorbance of each and every well ( $A/A0 \times 100$ ).







---

## 9. RECALLING AND RUNNING OF STORED TEST/PROGRAMS

The entire saved tests are available on the List Test or on the first screen after initialization. User can directly RUN the test by selecting TEST present on the screen.




After selection of Tests VIEW, DEL & RUN buttons are enable



### 9.1. Absorbance Mode:

Select the test and point to RUN option



On selecting RUN option, it will display following screen.



Again select run option present in run screen.

NAME ABSORB		MODE ABS		PRI 450		SEC 0						
SMPL 96		RUN										
	W1	W2	W3	W4	W5	W6	W7	W8	W9	W10	W11	W12
A	S1 0.000	S9 0.000	S17 0.000	S25 0.000	S33 0.000	S41 0.000	S49 0.000	S57 0.000	S65 0.000	S73 0.000	S81 0.000	S89 0.000
B	S2 0.000	S10 0.000	S18 0.000	S26 0.000	S34 0.000	S42 0.000	S50 0.000	S58 0.000	S66 0.000	S74 0.000	S82 0.000	S90 0.000
C	S3 0.000	S11 0.000	S19 0.000	S27 0.000	S35 0.000	S43 0.000	S51 0.000	S59 0.000	S67 0.000	S75 0.000	S83 0.000	S91 0.000
D	S4 0.000	S12 0.000	S20 0.000	S28 0.000	S36 0.000	S44 0.000	S52 0.000	S60 0.000	S68 0.000	S76 0.000	S84 0.000	S92 0.000
E	S5 0.000	S13 0.000	S21 0.000	S29 0.000	S37 0.000	S45 0.000	S53 0.000	S61 0.000	S69 0.000	S77 0.000	S85 0.000	S93 0.000
F	S6 0.000	S14 0.000	S22 0.000	S30 0.000	S38 0.000	S46 0.000	S54 0.000	S62 0.000	S70 0.000	S78 0.000	S86 0.000	S94 0.000
G	S7 0.000	S15 0.000	S23 0.000	S31 0.000	S39 0.000	S47 0.000	S55 0.000	S63 0.000	S71 0.000	S79 0.000	S87 0.000	S95 0.000
H	S8 0.000	S16 0.000	S24 0.000	S32 0.000	S40 0.000	S48 0.000	S56 0.000	S64 0.000	S72 0.000	S80 0.000	S88 0.000	S96 0.000
ESC		SAVE RESULT		LOAD NEXT		PRINT MATRX						

The instrument will ask “Plate Insert ? YES / NO”, select ‘Yes’. Instrument will read the absorbance using mechanical plate movement.

NAME ABSORB		MODE ABS		PRI 450		SEC 0						
SMPL 96		RUN										
	W1	W2	W3	W4	W5	W6	W7	W8	W9	W10	W11	W12
A	S1 0.000	S9 0.000	S17 0.000	S25 0.000	S33 0.000	S41 0.000	S49 0.000	S57 0.000	S65 0.000	S73 0.000	S81 0.000	S89 0.000
B	S2 0.000	S10 0.000	S18 0.000	S26 0.000	S34 0.000	S42 0.000	S50 0.000	S58 0.000	S66 0.000	S74 0.000	S82 0.000	S90 0.000
C	S3 0.000	S11 0.000	S19 0.000	S27 0.000	S35 0.000	S43 0.000	S51 0.000	S59 0.000	S67 0.000	S75 0.000	S83 0.000	S91 0.000
D	S4 0.000	S12 0.000	S20 0.000	S28 0.000	S36 0.000	S44 0.000	S52 0.000	S60 0.000	S68 0.000	S76 0.000	S84 0.000	S92 0.000
E	S5 0.000	S13 0.000	S21 0.000	S29 0.000	S37 0.000	S45 0.000	S53 0.000	S61 0.000	S69 0.000	S77 0.000	S85 0.000	S93 0.000
F	S6 0.000	S14 0.000	S22 0.000	S30 0.000	S38 0.000	S46 0.000	S54 0.000	S62 0.000	S70 0.000	S78 0.000	S86 0.000	S94 0.000
G	S7 0.000	S15 0.000	S23 0.000	S31 0.000	S39 0.000	S47 0.000	S55 0.000	S63 0.000	S71 0.000	S79 0.000	S87 0.000	S95 0.000
H	S8 0.000	S16 0.000	S24 0.000	S32 0.000	S40 0.000	S48 0.000	S56 0.000	S64 0.000	S72 0.000	S80 0.000	S88 0.000	S96 0.000
Plate Insert ?											YES	
											NO	

while reading the absorbance it will display the message “Wait for processing optical data....”.

NAME ABSORB		MODE ABS		PRI 450		SEC 0						
SMPL 96		RUN										
	W1	W2	W3	W4	W5	W6	W7	W8	W9	W10	W11	W12
A	S1 0.000	S9 0.000	S17 0.000	S25 0.000	S33 0.000	S41 0.000	S49 0.000	S57 0.000	S65 0.000	S73 0.000	S81 0.000	S89 0.000
B	S2 0.000	S10 0.000	S18 0.000	S26 0.000	S34 0.000	S42 0.000	S50 0.000	S58 0.000	S66 0.000	S74 0.000	S82 0.000	S90 0.000
C	S3 0.000	S11 0.000	S19 0.000	S27 0.000	S35 0.000	S43 0.000	S51 0.000	S59 0.000	S67 0.000	S75 0.000	S83 0.000	S91 0.000
D	S4 0.000	S12 0.000	S20 0.000	S28 0.000	S36 0.000	S44 0.000	S52 0.000	S60 0.000	S68 0.000	S76 0.000	S84 0.000	S92 0.000
E	S5 0.000	S13 0.000	S21 0.000	S29 0.000	S37 0.000	S45 0.000	S53 0.000	S61 0.000	S69 0.000	S77 0.000	S85 0.000	S93 0.000
F	S6 0.000	S14 0.000	S22 0.000	S30 0.000	S38 0.000	S46 0.000	S54 0.000	S62 0.000	S70 0.000	S78 0.000	S86 0.000	S94 0.000
G	S7 0.000	S15 0.000	S23 0.000	S31 0.000	S39 0.000	S47 0.000	S55 0.000	S63 0.000	S71 0.000	S79 0.000	S87 0.000	S95 0.000
H	S8 0.000	S16 0.000	S24 0.000	S32 0.000	S40 0.000	S48 0.000	S56 0.000	S64 0.000	S72 0.000	S80 0.000	S88 0.000	S96 0.000
Wait for processing optical data....												
ESC		SAVE RESULT		LOAD NEXT		PRINT MATRX						

After completion of reading it will ask: "Plate Remove? YES / NO". Remove the plate and Select 'YES'.

NAME ABSORB		MODE ABS		PRI 450		SEC 0						
SMPL 96		RUN										
	W1	W2	W3	W4	W5	W6	W7	W8	W9	W10	W11	W12
A	S1 0.000	S9 0.000	S17 0.000	S25 0.000	S33 0.000	S41 0.000	S49 0.000	S57 0.000	S65 0.000	S73 0.000	S81 0.000	S89 0.000
B	S2 0.000	S10 0.000	S18 0.000	S26 0.000	S34 0.000	S42 0.000	S50 0.000	S58 0.000	S66 0.000	S74 0.000	S82 0.000	S90 0.000
C	S3 0.000	S11 0.000	S19 0.000	S27 0.000	S35 0.000	S43 0.000	S51 0.000	S59 0.000	S67 0.000	S75 0.000	S83 0.000	S91 0.000
D	S4 0.000	S12 0.000	S20 0.000	S28 0.000	S36 0.000	S44 0.000	S52 0.000	S60 0.000	S68 0.000	S76 0.000	S84 0.000	S92 0.000
E	S5 0.000	S13 0.000	S21 0.000	S29 0.000	S37 0.000	S45 0.000	S53 0.000	S61 0.000	S69 0.000	S77 0.000	S85 0.000	S93 0.000
F	S6 0.000	S14 0.000	S22 0.000	S30 0.000	S38 0.000	S46 0.000	S54 0.000	S62 0.000	S70 0.000	S78 0.000	S86 0.000	S94 0.000
G	S7 0.000	S15 0.000	S23 0.000	S31 0.000	S39 0.000	S47 0.000	S55 0.000	S63 0.000	S71 0.000	S79 0.000	S87 0.000	S95 0.000
H	S8 0.000	S16 0.000	S24 0.000	S32 0.000	S40 0.000	S48 0.000	S56 0.000	S64 0.000	S72 0.000	S80 0.000	S88 0.000	S96 0.000
Plate Remove?											YES	
											NO	

Details of absorbance of all wells for each strip will be displayed on the screen as shown in below

NAME ABSORB		MODE ABS		PRI 450		SEC 0						
SMPL 96		RUN										
	W1	W2	W3	W4	W5	W6	W7	W8	W9	W10	W11	W12
A	S1 0.028	S9 0.145	S17 0.100	S25 0.017	S33 0.089	S41 0.135	S49 0.362	S57 0.078	S65 0.410	S73 1.440	S81 0.746	S89 0.224
B	S2 0.058	S10 0.151	S18 0.386	S26 2.553	S34 0.111	S42 1.479	S50 2.593	S58 0.081	S66 2.517	S74 0.781	S82 0.730	S90 0.005
C	S3 0.065	S11 2.596	S19 0.135	S27 0.142	S35 0.141	S43 0.182	S51 0.108	S59 0.112	S67 1.585	S75 1.530	S83 0.405	S91 2.653
D	S4 0.072	S12 0.150	S20 2.630	S28 0.143	S36 2.650	S44 0.077	S52 0.159	S60 0.170	S68 1.504	S76 0.071	S84 0.426	S92 0.222
E	S5 2.541	S13 0.140	S21 0.149	S29 0.091	S37 0.166	S45 0.173	S53 2.653	S61 0.191	S69 0.097	S77 2.659	S85 0.252	S93 0.222
F	S6 2.614	S14 2.682	S22 0.185	S30 0.130	S38 2.648	S46 0.167	S54 0.176	S62 0.071	S70 2.686	S78 0.853	S86 0.907	S94 0.426
G	S7 0.056	S15 0.163	S23 2.571	S31 0.155	S39 0.147	S47 0.119	S55 2.618	S63 0.081	S71 0.085	S79 0.070	S87 0.407	S95 2.591
H	S8 0.071	S16 0.123	S24 0.113	S32 0.054	S40 0.116	S48 0.060	S56 0.112	S64 0.111	S72 0.093	S80 2.568	S88 0.396	S96 0.015
ESC		SAVE RESULT		LOAD NEXT		PRINT MATRX						

In this operation if you keep BLANK YES.

The first well of first strip will be considered as blank and instrument will give the absorbance of remaining all well with blank subtraction from original absorbance.

The “PRINT MATRIX” will print the details of entire plate with well identification and absorbance of that well, like.....

A

W1    W2    W3    W4    W5    W6  
0.125  0.250  0.350  0.450  0.550  0.650..

“SAVE RESULT” is used to save the result.

“LOAD NEXT” is used to load the Next Plate in order to get the absorbance of next plate

**Example Print-out of ABSORBANCE test obtained from micro read 1000**

After running the plate in Absorbance mode, user can take the print of the result in the form of print matrix which will be printed as follows –

```

ABS,          21/07/08          11:47:34
A
W1   W2   W3   W4   W5   W6   W7   W8   W9   W10
0.028 0.145 0.100 0.017 0.089 0.135 0.362 0.078 0.410 1.441
B
W1   W2   W3   W4   W5   W6   W7   W8   W9   W10
0.058 0.151 0.386 2.553 0.111 1.479 2.593 0.081 2.517 0.781
C
W1   W2   W3   W4   W5   W6   W7   W8   W9   W10
0.065 2.596 0.135 0.142 0.141 0.182 0.108 0.112 1.585 1.531
D
W1   W2   W3   W4   W5   W6   W7   W8   W9   W10
0.072 0.150 2.630 0.143 2.650 0.077 0.159 0.170 1.504 0.071
E
W1   W2   W3   W4   W5   W6   W7   W8   W9   W10
2.541 0.140 0.149 0.091 0.166 0.173 2.653 0.191 0.097 2.651
F
W1   W2   W3   W4   W5   W6   W7   W8   W9   W10
2.614 2.682 0.185 0.130 2.648 0.167 0.176 0.071 2.686 0.851
G
W1   W2   W3   W4   W5   W6   W7   W8   W9   W10
0.058 0.145 0.100 0.017 0.089 0.135 0.362 0.078 0.410 1.441

```



## 9.2. Cut off Mode

User can be recalling the saved test in Cut-Off mode by selecting view, it will display following screen.

User can take the print of the test parameters by selecting “PRN” option present on the screen. It will print all the details of the test parameters along with QC equations.

NAME TSH	MODE COFF	PRI 405	SEC 630	ESC
BLK Y 1 0.100	NC Y 3 0.150	PC Y 2 0.000	LC N 0 0.000	SAVE
CC N 0 0.000	QC 0			
CUT-OFF INDEX 0.000 0.900		GREY ZONE 0.000 0.000		
CUT-OFF EQUATION NC+PC				

A) For running the test in CUTOFF mode select the particular test and point to RUN option present on the “List Test” screen.

SIN							
TSH							
UPT							
CPC							
KIN							
CRP							
NXT	PRV	NEW	VIEW	DEL	MENU	RUN	PRN



S

NAME TSH		MODE COFF		PRI 450		SEC 0		BLANK 1		NC 3		PC 2	
SMPL 0		RUN		PID		GREYZONE 10.00				CABS			
	W1	W2	W3	W4	W5	W6	W7	W8	W9	W10	W11	W12	
A	B1 0.000												
B	NC1 0.000												
C	NC2 0.000												
D	NC3 0.000												
E	PC1 0.000												
F	PC2 0.000												
G													
H													
										SIMUL		MODIFY	
ESC		SAVE RESULT		LOAD NEXT		PRINT MATRX		PRINT RESULT		ACCEPT TEST			

**B) It shows the 8 strips vertical with 12 wells Horizontal with blanks and controls loaded in first strip. Select the ‘SMPL and enter the total Number of samples you want to load with the help of numeric screen displayed at the bottom of the screen.**

NAME TSH		MODE COFF		PRI 450		SEC 0		BLANK 1		NC 3		PC 2					
SMPL 0		RUN		PID		GREYZONE 10.00				CUT-OFF ABSORBANCE							
	W1	W2	W3	W4	W5	W6	W7	W8	W9	W10	W11	W12					
A	B1 0.000																
B	NC1 0.000																
C	NC2 0.000																
D	NC3 0.000																
E	PC1 0.000																
F	PC2 0.000																
G																	
H																	
SMPL 0		DUP N		ENTER NUM:													
ESC		SAVE RESULT		1		2		3		4		5		6		ENT	
				7		8		9		0		.		CLR			

After entering number of sample the run screen will be display in the below

NAME TSH		MODE COFF		PRI 450		SEC 0		BLANK 1		NC 3		PC 2	
SMPL 10		RUN		PID		GREYZONE 10.00				CUT-OFF ABSORBANCE			
	W1	W2	W3	W4	W5	W6	W7	W8	W9	W10	W11	W12	
A	B1 0.000	S3 0.000											
B	NC1 0.000	S4 0.000											
C	NC2 0.000	S5 0.000											
D	NC3 0.000	S6 0.000											
E	PC1 0.000	S7 0.000											
F	PC2 0.000	S8 0.000											
G	S1 0.000	S9 0.000											
H	S2 0.000	S10 0.000											
										SIMUL		MODIFY	
ESC		SAVE RESULT		LOAD NEXT		PRINT MATRX		PRINT RESULT		ACCEPT TEST			

C) After entering No. of samples, select RUN. The plate-loading tray will come out and the following message will appear: “Plate Insert? YES / NO”

D) Here the plate is loaded in the direction that controls should go in first. Select “YES”.

NAME TSH		MODE COFF		PRI 450		SEC 0		BLANK 1		NC 3		PC 2	
SMPL 0		RUN		PID		GREYZONE 10.00				CUT-OFF ABSORBANCE			
	W1	W2	W3	W4	W5	W6	W7	W8	W9	W10	W11	W12	
A	B1 0.000	S3 0.000											
B	NC1 0.000	S4 0.000											
C	NC2 0.000	S5 0.000											
D	NC3 0.000	S6 0.000											
E	PC1 0.000	S7 0.000											
F	PC2 0.000	S8 0.000											
G	S1 0.000	S9 0.000											
H	S2 0.000	S10 0.000											
Plate Insert ?												YES	
												NO	
ESC		SAVE RESULT		LOAD NEXT		PRINT MATRX		PRINT RESULT		ACCEPT TEST			

The instrument will read the controls and samples while displaying the message “Wait for processing optical data...”.

NAME TSH		MODE COFF		PRI 450		SEC 0		BLANK 1		NC 3		PC 2	
SMPL 0		RUN		PID		GREYZONE 10.00				CUT-OFF ABSORBANCE			
W1	W2	W3	W4	W5	W6	W7	W8	W9	W10	W11	W12		
A	B1 0.000	S3 0.000											
B	NC1 0.000	S4 0.000											
C	NC2 0.000	S5 0.000											
D	NC3 0.000	S6 0.000											
E	PC1 0.000	S7 0.000											
F	PC2 0.000	S8 0.000											
G	S1 0.000	S9 0.000											
H	S2 0.000	S10 0.000											

SIMUL      MODIFY

Wait for processing optical data...

ESC	SAVE RESULT	LOAD NEXT	PRINT MATRX	PRINT RESULT	ACCEPT TEST
-----	----------------	--------------	----------------	-----------------	----------------

**Following message appears: “Plate Remove? YES / NO” after optical data processing completed.  
Remove the plate and select YES.**

NAME TSH		MODE COFF	PRI 450	SEC 0	BLANK 1	NC 3	PC 2				
SMPL 0	RUN	PID	GREYZONE 10.00			CUT-OFF ABSORBANCE					
W1	W2	W3	W4	W5	W6	W7	W8	W9	W10	W11	W12
A	B1 0.000	S3 0.000									
B	NC1 0.000	S4 0.000									
C	NC2 0.000	S5 0.000									
D	NC3 0.000	S6 0.000									
E	PC1 0.000	S7 0.000									
F	PC2 0.000	S8 0.000									
G	S1 0.000	S9 0.000									
H	S2 0.000	S10 0.000									
Plate Remove ?										YES	
										NO	
ESC	SAVE RESULT	LOAD NEXT	PRINT MATRX	PRINT RESULT	ACCEPT TEST						

E) After selecting YES, result screen will appear with the measured values of blank, Controls and the cut off absorbance.

NAME TSH		MODE COFF		PRI 450		SEC 0		BLANK 1		NC 3		PC 2	
SMPL 0		RUN		PID		GREYZONE 10.00				CUT-OFF ABSORBANCE 0.180			
	W1	W2	W3	W4	W5	W6	W7	W8	W9	W10	W11	W12	
A	B1 0.002	S3 0.145											
B	NC1 0.082	S4 0.000											
C	NC2 0.040	S5 2.573											
D	NC3 0.080	S6 0.157											
E	PC1 0.080	S7 0.172											
F	PC2 2.682	S8 2.612											
G	S1 2.626	S9 0.091											
H	S2 0.147	S10 0.382											
								TABLE		SIMUL		MODIFY	
ESC		SAVE RESULT		LOAD NEXT		PRINT MATRX		PRINT RESULT		ACCEPT TEST			



---

**G)** By selecting ‘ACCEPT TEST’ option you can save this test with the details of controls absorbance. Next time it is possible to run the same test without loading controls. This means that the previous stored values of controls can be used.

All the details are in table format, as per row and column shown on the screen. You can obtain the print of the same format by selecting option ‘PRINT RESULT’.

You can obtain the print in matrix format, by selecting option ‘PRINT MATRX’.

In matrix form you will get the print as per your tray/plate for all the wells. Instrument prints seven different parameters in print matrix for a single well.

ROW IDENTIFICATION	A, B, C,D,.....
WELL NUMBER	W1, W2, W3, W4.....
WELL ID OR SAMPLE NUMBER	BL, NC, PC, OR S1, S2, S3, .....
PATIENT IDENTIFICATION	RAMESH, VIKAS, etc.
WELL ABSORBANCE	0.050, 0.098, 0.085, ...
SAMPLE CONCENTRATION	0.738, 0.689, 2.578,...
INTERPRETATION / REMARKS	POS or NEG or EQ

**Example Print-out of POS-CUTOFF test obtained from micro read 1000**

TSH, 21/07/08,16:03:08  
 Range:- POS>0.087, NEG<0.071

Cutabs=0.0793

<b>A</b>											
W1	W2	W3	W4	W5	W6	W7	W8	W9	W10	W11	W12
B	S4	S12	S20	S28	S36	S44	S52	S60	S68	S76	S84
	OPD13	OPD21	OPD29	OPD37	OPD45	OPD53	OPD61	OPD69	OPD77	OPD85	OPD93
0.002	0.147	0.091	0.011	0.083	0.127	0.362	0.072	0.400	1.454	0.743	0.222
0.000	1.302	0.803	0.099	0.737	1.122	3.201	0.836	3.533	12.84	6.564	1.959
	POS	NEG	NEG	NEG	POS	POS	NEG	POS	POS	POS	POS
<b>B</b>											
W1	W2	W3	W4	W5	W6	W7	W8	W9	W10	W11	W12
NC	S5	S13	S21	S29	S37	S45	S53	S61	S69	S77	S85
	OPD14	OPD22	OPD30	OPD38	OPD46	OPD54	OPD62	OPD70	OPD78	OPD86	OPD94
0.062	0.145	0.382	2.535	0.116	1.470	2.573	0.115	2.515	0.747	0.748	0.003
0.000	1.278	3.372	22.39	1.024	12.99	22.73	1.011	22.22	6.597	6.608	0.025
	POS	POS	POS	EQ	POS	POS	EQ	POS	POS	POS	NEG
<b>C</b>											
W1	W2	W3	W4	W5	W6	W7	W8	W9	W10	W11	W12
NC	S6	S14	S22	S30	S38	S46	S54	S62	S70	S78	S86
	OPD15	OPD23	OPD31	OPD39	OPD47	OPD55	OPD63	OPD71	OPD79	OPD87	OPD95
0.040	2.573	0.129	0.159	0.133	0.217	0.077	0.091	1.584	1.557	0.502	2.638
0.000	22.72	1.136	1.405	1.177	1.917	0.677	0.804	13.99	13.75	4.434	23.30
	POS	POS	POS	POS	POS	NEG	NEG	POS	POS	POS	POS
<b>D</b>											
W1	W2	W3	W4	W5	W6	W7	W8	W9	W10	W11	W12
PC	S7	S15	S23	S31	S39	S47	S55	S63	S71	S79	S87
	OPD16	OPD24	OPD32	OPD40	OPD48	OPD56	OPD64	OPD72	OPD80	OPD88	OPD96
0.080	0.157	2.630	0.141	2.642	0.058	0.155	0.159	1.515	0.071	0.423	0.209
0.000	1.389	23.23	1.246	23.33	0.516	1.370	1.403	13.38	0.628	3.732	1.850
	POS	POS	POS	POS	NEG	POS	POS	POS	NEG	POS	POS
<b>E</b>											
W1	W2	W3	W4	W5	W6	W7	W8	W9	W10	W11	W12
PC	S8	S16	S24	S32	S40	S48	S56	S64	S72	S80	S88
	OPD17	OPD25	OPD33	OPD41	OPD49	OPD57	OPD65	OPD73	OPD81	OPD89	OPD97
2.562	0.172	0.198	0.061	0.182	0.165	2.656	0.156	0.080	2.616	0.259	0.224
0.000	1.517	1.751	0.536	1.612	1.457	23.46	1.377	0.705	23.11	2.291	1.978
	POS	POS	NEG	POS	POS	POS	POS	NEG	POS	POS	POS
<b>F</b>											
W1	W2	W3	W4	W5	W6	W7	W8	W9	W10	W11	W12
S1	S9	S17	S25	S33	S41	S49	S57	S65	S73	S81	S89
OPD10	OPD18	OPD26	OPD34	OPD42	OPD50	OPD58	OPD66	OPD74	OPD82	OPD90	OPD98
2.626	2.612	0.055	0.116	2.649	0.188	0.176	0.052	2.638	0.855	0.899	0.419
23.19	23.07	0.483	1.023	23.40	1.660	1.553	0.460	23.30	7.548	7.942	3.699
	POS	NEG	EQ	POS	POS	POS	NEG	POS	POS	POS	POS
<b>G</b>											
W1	W2	W3	W4	W5	W6	W7	W8	W9	W10	W11	W12
S2	S10	S18	S26	S34	S42	S50	S58	S66	S74	S82	S90
OPD11	OPD19	OPD27	OPD35	OPD43	OPD51	OPD59	OPD67	OPD75	OPD83	OPD91	OPD99
0.049	0.140	2.595	0.138	0.142	0.117	2.611	0.076	0.081	0.050	0.403	2.579
0.429	1.234	22.92	1.215	1.252	1.035	23.06	0.672	0.717	0.441	3.564	22.78
	NEG	POS	POS	POS	EQ	POS	NEG	NEG	NEG	POS	POS
<b>H</b>											
W1	W2	W3	W4	W5	W6	W7	W8	W9	W10	W11	W12
S3	S11	S19	S27	S35	S43	S51	S59	S67	S75	S83	S91
OPD12	OPD20	OPD28	OPD36	OPD44	OPD52	OPD60	OPD68	OPD76	OPD84	OPD92	OPD100
0.061	0.109	0.002	0.071	0.076	0.131	0.109	0.101	0.082	2.560	0.399	0.015
0.538	0.964	0.018	0.631	0.670	1.157	0.960	1.896	0.721	22.61	3.527	0.132
	NEG	EQ	NEG	NEG	POS	EQ	NEG	NEG	POS	POS	NEG

### 9.2. A. Invalid Assay in Cut Off Mode:

If any individual control behaves incorrectly and its absorbance is not satisfying the QC check value of that control, the following message appears: “Invalid Assay” and a remark “HI” or “LO” for that particular control. OR if the QC1, QC2, QC3 and QC4, which has an other QC condition of controls (like some reagent manual gives condition  $PC-NC > 0.2$  or  $NC/PC > 0.5$ ) the following message appears: “Invalid Assay” and a remark either ‘HI’ or ‘LO’ for that QC condition.

If “QC” is selected, then first screen is display after running the test is shown in below

NAME TSH		PRI 450	SEC 450	CABS 0.118	R POS 0.000	R NEG 0.000	
BLK 1	QC Y	VAL 0.20	AVG 0.002	REM	QC1 0.000	REM	ESC
NC 3	QC Y	VAL 0.25	AVG 0.233	REM HI	QC2 0.000	REM	PRN
PC 2	QC Y	VAL 0.35	AVG 0.150	REM	QC3 0.000	REM	
LC 0	QC N	VAL	AVG	REM	QC4 0.000	REM	
CC 0	QC N	VAL	AVG	REM			

Invalid Assay

If "QC" is not selected

NAME TSH		MODE COFF		PRI 450		SEC 0		BLANK 1		NC 3		PC 2	
SMPL 0		RUN		PID		GREYZONE 10.00				CUT-OFF ABSORBANCE 0.180			
	W1	W2	W3	W4	W5	W6	W7	W8	W9	W10	W11	W12	
A	B1 0.002	S3 3.600											
B	NC1 0.082	S4 3.670											
C	NC2 0.040	S5 2.573											
D	NC3 0.080	S6 0.157											
E	PC1 0.080	S7 0.172											
F	PC2 2.682	S8 2.612											
G	S1 2.626	S9 0.091											
H	S2 0.147	S10 0.382											
								TABLE		SIMUL		MODIFY	
ESC		SAVE RESULT		LOAD NEXT		PRINT MATRX		PRINT RESULT		ACCEPT TEST			


(NOTE: The absorbance of controls should not go above 3.50. If it is > than 3.50, it will show the remark 'HI' for that control. In such a case edit the control absorbance.)

In such case, delete an individual control, so that the average of remaining controls will satisfy the QC condition of that control as per reagent manual.

To select any control for deleting, you must select the touch zone as shown in the below screen.

When you select any particular control for deleting (suppose you have selected Negative Control (NC), 'AVG' appear on the screen.

“Select DEL key to proceed”. To delete the control, select 'DEL'. After deleting that control, the absorbance of that control will become zero and the average of controls will also get changed. Select 'ESC' option to come back to the QC Result screen.



When you go back, you will get the edited values for average of control with no remark for any control and no “Invalid Assay” message.

**\*NOTE:**

Make the blank absorbance valid in case the absorbance doesn't meet the QC conditions.

If necessary valid the control absorbance because the absorbance of the blank affect the absorbance of all controls and samples.

It is recommended that, for a test with single blank and single control, if absorbance doesn't satisfy QC conditions or if all controls of any single type (means all NC controls or all PC controls) doesn't meet the QC, the test will become totally invalid. It is not possible to make the VALID TEST and new controls have to be loaded.

### 9.3. Multistandard:

Select the multistandard mode test from the Test List screen.

#### Multi Standard Mode

NAME MUL_4P	MODE MSTD	PRI 450	SEC 0	ESC
BLK N Y 0.000	CAL DUP 7	HI CTRL 1 25.03 10.54	LO CTRL 1 16.61 5.45	EDIT
NORMAL RANGE 15.25 9.46		UNIT mg/ml		PRN

A) For running the test in MSTD mode select the particular test and point to RUN option present on the Test list screen.

SIN							
TSH							
UPT							
CPC							
KIN							
CRP							
NXT	PRV	NEW	VIEW	DEL	MENU	RUN	PRN

On selecting RUN option it will display following screen.

NAME TSH		MODE MSTD		PRI 450		SEC 0		BLK Y		CALI 7		HC 1		LC 1	
SMPL 2		RUN		PID		HI CTRL 10.540-25.030			LOW CTRL 5.450-16.610			RANGE 9.460-15.250			
	W1	W2	W3	W4	W5	W6	W7	W8	W9	W10	W11	W12			
A	B1 0.000														
B	C1 0.000														
C	C2 0.000														
D	C3 0.000														
E	C4 0.000														
F	C5 0.000														
G	C6 0.000														
H	C7 0.000														
ESC		SAVE RESULT		LOA NEX		NXT	PRV	NEW	VIEW	DEL	MENU	RUN	PRN		

SIN

TSH

UPT

CPC

KIN

CRP

Select the empty well where you have to locate high control (HC1) and then select “HICO” which shown in following screen.

NAME TSH		MODE MSTD		PRI 450		SEC 0		BLK Y		CALI 7		HC 1		LC 1	
SMPL 2		RUN		PID		HI CTRL 10.540-25.030				LOW CTRL 5.450-16.610			RANGE 9.460-15.250		
	W1	W2	W3	W4	W5	W6	W7	W8	W9	W10	W11	W12			
A	B1 0.000														
B	C1 0.000														
C	C2 0.000														
D	C3 0.000														
E	C4 0.000														
F	C5 0.000														
G	C6 0.000														
H	C7 0.000														

HICO	LOCO	ESC	SIMUL	MODIFY
------	------	-----	-------	--------

RESULT	NEXT	MATRX	RESULT	ACCEPT TEST	GRAPH
--------	------	-------	--------	-------------	-------



After loading high control, screen is display below.

NAME TSH		MODE MSTD	PRI 450	SEC 0	BLK Y	CALI 7	HC 1	LC 1				
SMPL 2	RUN	PID	HI CTRL 10.540-25.030			LOW CTRL 5.450-16.610		RANGE 9.460-15.250				
W1	W2	W3	W4	W5	W6	W7	W8	W9	W10	W11	W12	
A	B1 0.000	HC1 0.000										
B	C1 0.000											
C	C2 0.000											
D	C3 0.000											
E	C4 0.000											
F	C5 0.000											
G	C6 0.000											
H	C7 0.000											
										SIMUL	MODIFY	
ESC	SAVE RESULT	LOAD NEXT	PRINT MATRX	PRINT RESULT	ACCEPT TEST	GRAPH						

Similarly, you can enter No. of LO Controls and locate its positions.

User can change the position of blank, calibrator & controls(High/low) with help of “MODIFY”

Which shown in following fig

NAME TSH		MODE MSTD	PRI 450	SEC 0	BLK Y	CALI 7	HC 1	LC 1				
SMPL 2	RUN	PID	HI CTRL 10.540-25.030			LOW CTRL 5.450-16.610		RANGE 9.460-15.250				
W1	W2	W3	W4	W5	W6	W7	W8	W9	W10	W11	W12	
A	B1 0.000	HC1 0.000										
B	C1 0.000	LC1 0.000										
C	C2 0.000											
D	C3 0.000											
E	C4 0.000											
F	C5 0.000											
G	C6 0.000											
H	C7 0.000											
										SIMUL	MODIFY	
ESC	SAVE RESULT	LOAD NEXT	PRINT MATRX	PRINT RESULT	ACCEPT TEST	GRAPH						

After selecting “MODIFY”, pop up message is display i.e. “Select Well to load blank/ Calibrator”

Which shown in following fig

NAME TSH		MODE MSTD		PRI 450		SEC 0		BLK Y		CALI 7		HC 1		LC 1	
SMPL 2		RUN		PID		HI CTRL 10.540-25.030			LOW CTRL 5.450-16.610			RANGE 9.460-15.250			
W1	W2	W3	W4	W5	W6	W7	W8	W9	W10	W11	W12				
A															
B															
C															
D															
E															
F															
G															
H															
											SIMUL		MODIFY		
Select well to load blank/ calibrator															
ESC		SAVE RESULT		LOAD NEXT		PRINT MATRX		PRINT RESULT		ACCEPT TEST		GRAPH			

After pop up message, select the well where you want to load the blank or calibrator or controls

**For e.g. Suppose you select column 2 & row A for HI control will be display on screen as shown below.**

NAME TSH		MODE MSTD		PRI 450		SEC 0		BLK Y		CALI 7		HC 1		LC 1	
SMPL 2		RUN		PID		HI CTRL 10.540-25.030				LOW CTRL 5.450-16.610		RANGE 9.460-15.250			
	W1	W2	W3	W4	W5	W6	W7	W8	W9	W10	W11	W12			
A															
B															
C															
D															
E															
F															
G															
H															
										SIMUL		MODIFY			
ESC		SAVE RESULT		LOAD NEXT		PRINT MATRX		PRINT RESULT		ACCEPT TEST		GRAPH			

NAME TSH		MODE MSTD		PRI 450		SEC 0		BLK Y		CALI 7		HC 1		LC 1	
SMPL 2		RUN		PID		HI CTRL 10.540-25.030				LOW CTRL 5.450-16.610		RANGE 9.460-15.250			
W1	W2	W3	W4	W5	W6	W7	W8	W9	W10	W11	W12				
A	HC1 0.000														
B															
C															
D															
E															
F															
G															
H															
										SIMUL		MODIFY			
ESC		SAVE RESULT		LOAD NEXT		PRINT MATRX		PRINT RESULT		ACCEPT TEST		GRAPH			

Repeat the above step to load remaining calibrators, control and blank.

After loading blank, calibrators & low controls screen will be display shown in below.

NAME TSH		MODE MSTD	PRI 450	SEC 0	BLK Y	CALI 7	HC 1	LC 1				
SMPL 2	RUN	PID	HI CTRL 10.540-25.030			LOW CTRL 5.450-16.610		RANGE 9.460-15.250				
	W1	W2	W3	W4	W5	W6	W7	W8	W9	W10	W11	W12
A	B1 0.000	HC1 0.000										
B	C1 0.000	LC1 0.000										
C	C2 0.000											
D	C3 0.000											
E	C4 0.000											
F	C5 0.000											
G	C6 0.000											
H	C7 0.000											

SIMUL MODIFY

ESC SAVE RESULT LOAD NEXT PRINT MATRX PRINT RESULT ACCEPT TEST GRAPH

E) Select “SMPL” you want to run. In case of 10 samples, it will load 10 samples as follows -

NAME TSH		MODE MSTD	PRI 450	SEC 0	BLK Y	CALI 7	HC 1	LC 1			
SMPL 2	RUN	PID	HI CTRL 10.540-25.030			LOW CTRL 5.450-16.610		RANGE 9.460-15.250			
W1	W2	W3	W4	W5	W6	W7	W8	W9	W10	W11	W12
A	B1 0.000	HC1 0.000	S7 0.000								
B	C1 0.000	LC1 0.000	S8 0.000								
C	C2 0.000	S1 0.000	S9 0.000								
D	C3 0.000	S2 0.000	S10 0.000								
E	C4 0.000	S3 0.000									
F	C5 0.000	S4 0.000									
G	C6 0.000	S5 0.000									
H	C7 0.000	S6 0.000									

**F) Now Select RUN option so that readings can be taken in MULTI STANDARD mode.**

The plate holder will come out and will display the message string- “Plate Insert? YES / NO”. The plate has to be loaded in such a way that the direction of the calibrators should go in first and then select ‘YES’.

The instrument will read the absorbance of calibrators and samples by displaying the message string “wait for processing optical data...”. And finally display the string- “Plate Remove? YES / NO”. Remove the plate and select ‘YES’. The plate holder will go inside and following screen will be flashed on the display.

NAME		MODE	PRI	SEC	BLK	CALI	HC	LC				
TSH		MSTD	450	0	Y	7	1	1				
SMPL	RUN	PID	HI CTRL			LOW CTRL			RANGE			
2			10.540-25.030			5.450-16.610			9.460-15.250			
W1	W2	W3	W4	W5	W6	W7	W8	W9	W10	W11	W12	
A	B1 0.050	HC1 0.350	S7 0.222									
B	C1 0.070	LC1 0.248	S8 0.500									
C	C2 0.257	S1 0.200	S9 0.920									
D	C3 0.280	S2 1.854	S10 2.050									
E	C4 0.285	S3 1.005										
F	C5 0.295	S4 1.250										
G	C6 0.300	S5 1.500										
H	C7 0.320	S6 1.925										
						TABLE	SIMUL	MODIFY				
ESC	SAVE RESULT	LOAD NEXT	PRINT MATRX	PRINT RESULT	ACCEPT TEST	GRAPH						

---

All the details as per row and column are shown on the screen. The printout of the same format can be obtained by selecting 'PRINT RESULTS'.

You can get the print in matrix format, by selecting 'PRINT MATRX'.

In matrix form you will get the print as per your plate for all the wells. Instrument will print seven different parameters in print matrix for a single well as follows:

ROW IDENTIFICATION	A,B, C,D,.....
WELL NUMBER	W1, W2, W3, W4...
WELL ID OR SAMPLE NUMBER	B, C1, C2, C3, OR S1, S2, S3,..
PATIENT IDENTIFICATION	RAMESH, VIKAS, etc.
WELL ABSORBANCE	0.050, 0.098, 0.085, ...
SAMPLE CONCENTRATION	0.738, 0.689, 2.578...
INTERPRETATION / REMARKS	POS, NEG, EQ,

**Remember that Negative sign result is interpreted as 0.001 both in printout as well as on screen.**




**Example Print-out of MultiStandard test obtained from micro read 1000**

```

MULTI-S
Test Name :- MUL
Pri:450 Sec:000
BL N          DUP N
CAL 05        DUP N
UNIT:ug/ml
HI CO:High >N      Low
LO CO:High >N      Low
Range:High >6.000  Low
Graph: PT TO PT
Y vs X: ABS VS CON
No      Con      Abs
01      10.000    0.079
02      17.000    0.084
    
```

X=0.625, Y=0.126




MUL. 29/12/19, 12:59:52,  
Range :- HI > 6:000, LO < 4:000

A											
U1	U2	U3	U4	U5	U6	U7	U8	U9	U10	U11	U12
C1	S4	S12	S20	S28	S36	S44	S52	S60	S68	S76	S84
	4	12	20	28	36	44	52	60	68	76	84
1.259	1.218	0.748	1.499	1.324	1.885	1.629	1.284	2.001	1.566	4.199	2.536
00010	00013	00399	-7.32	05.29	-35.3	-16.7	000.2	-67.5	-12.2	-0200	-82.0
	ug/ml	ug/ml	ug/ml	ug/ml	ug/ml	ug/ml	ug/ml	ug/ml	ug/ml	ug/ml	ug/ml
	HI	HI	LO	EQU	LO	LO	HI	LO	LO	LO	LO

B											
U1	U2	U3	U4	U5	U6	U7	U8	U9	U10	U11	U12
C2	S5	S13	S21	S29	S37	S45	S53	S61	S69	S77	S85
	5	13	21	29	37	45	53	61	69	77	85
1.163	1.229	0.731	1.540	1.318	1.839	1.560	1.290	2.293	1.560	3.113	2.377
00017	012.2	00418	-18.9	05.72	-31.9	-11.8	07.78	-64.7	-11.8	-0124	-78.8
	ug/ml	ug/ml	ug/ml	ug/ml	ug/ml	ug/ml	ug/ml	ug/ml	ug/ml	ug/ml	ug/ml
	HI	HI	LO	EQU	LO	LO	HI	LO	LO	LO	LO

**G) On selecting 'Graph:', it displays the Graph screen with all the details of X-axis and Y-axis. Select 'PRINT' option to obtain the print out.**



Select 'ACCEPT TEST' option to save this test with the details of calibrator absorbance. It is possible to run the same test without loading calibrators. The previous graph can be used for new samples.

(\* **NOTE:** If you do not want to run the calibrators each time, please select "ACCEPT TEST" option to save graph, whenever you run the test with calibrators.)

### 9.3.A.Invalid assay in Multi standard Mode



In Multi standard, the absorbance of calibrators should be in increasing or decreasing from one calibrator to the next calibrator. If any calibrator/s behave(s) incorrectly you will obtain a message “Invalid Assay” at the bottom of the screen.

NAME TSH		MODE MSTD	PRI 450	SEC 0	BLK Y	CALI 7	HC 1	LC 1			
SMPL 2	RUN	PID	HI CTRL 10.540-25.030			LOW CTRL 5.450-16.610		RANGE 9.460-15.250			
W1	W2	W3	W4	W5	W6	W7	W8	W9	W10	W11	W12
A	B1 0.050	HC1 0.350	S7 0.222								
B	C1 0.070	LC1 0.248	S8 0.500								
C	C2 0.257	S1 0.200	S9 0.920								
D	C3 0.280	S2 1.854	S10 2.050								
E	C4 0.285	S3 1.005									
F	C5 0.295	S4 1.250									
G	C6 0.300	S5 1.500									
H	C7 0.320	S6 1.925									

TABLE	SIMUL	MODIFY
-------	-------	--------

ESC	SAVE RESULT	LOAD NEXT	PRINT MATRX	PRINT RESULT	ACCEPT TEST	GRAPH
-----	----------------	--------------	----------------	-----------------	----------------	-------

If Assay is Invalid then “ACCEPT TEST”, “SAVE RESULTS”, “PRINT MATRX”, “SAVE RESULT” options are blocked. In such a case, select GRAPH option and change the ASSAY to VALID with the help of EDIT option present on the Graph screen.



You can easily select any particular calibrator by touching the particular touch zone and edit the absorbance of that. After editing the calibrator’s absorbance select “ESC”. Screen will be blink and then the modified graph will appear on the screen.

If calibrator's value being edited in Multistandard mode then "\*Modified..." string is displayed in the print and in the run screen.

## 10. RERUNNING ACCEPTED TESTS / PROGRAMS

“ACCEPT TEST” option is used to store the data of controls or calibrators. After running any test first time with controls or calibrators, user can select the option “Accept Test” to save the data of controls or calibrators, so that next time when you want to load the same test, there is no need to load controls or calibrators in ELISA plate. You can use previously stored data. If user recall the test with saved parameters then "Using stored values" strings should appeared in the printout.

### A) In “Multi standard” Mode:

After Run the test select Accept test option present on run screen.

NAME TSH		MODE MSTD		PRI 450		SEC 0		BLK Y		CALI 7			
SMPL 2		RUN		PID		RANGE 10.30-16.81							
	W1	W2	W3	W4	W5	W6	W7	W8	W9	W10	W11	W12	
A	B1 0.050												
B	C1 0.070												
C	C2 0.257												
D	S1 0.200												
E	S2 1.854												
F													
G													
H													
								TABLE		SIMUL		MODIFY	
ESC		SAVE RESULT		LOAD NEXT		PRINT MATRX		PRINT RESULT		ACCEPT TEST		GRAPH	

After selecting Accept test option, it will display message i.e. “Test Accepted” .

The screenshot displays the F10 software interface with the following elements:



- Parameters:** NAME TSH, MODE MSTD, PRI 450, SEC 0, BLK Y, CALI 7.
- Operational Data:** SMPL 2, RUN, PID, RANGE 10.30-16.81.
- Table:** A grid with 12 columns (W1-W12) and 8 rows (A-H). The first two rows (A and B) contain data for B1 (0.050) and C1 (0.070) respectively. Rows C through H are currently empty.
- Message:** "Test Accepted..."
- Buttons:** TABLE, SIMUL, MODIFY, ESC, SAVE RESULT, LOAD NEXT, PRINT MATRX, PRINT RESULT, ACCEPT TEST, GRAPH.

	W1	W2	W3	W4	W5	W6	W7	W8	W9	W10	W11	W12
A	B1 0.050											
B	C1 0.070											
C	C2 0.257											
D	S1 0.200											
E	S2 1.854											
F												
G												
H												

**B) In Multi Standard / Single Standard / % Absorbance / Uptake/Ratio Mode:**

When you want to run such an accepted test, select that particular test for instant 'MULTISTANDARD' mode and point to 'RUN' option present on the Test screen.

following screen display.



S  
 ↓

NAME TSH		MODE MSTD		PRI 450		SEC 0		BLK Y		CALI 2		
SMPL 2		RUN		PID		RANGE 10.30-16.61						
	W1	W2	W3	W4	W5	W6	W7	W8	W9	W10	W11	W12
A												
B												
C												
D												
E												
F												
G												
H												
LOAD NEW CALIBRATOR?										YES		
										NO		

**NOTE:** User can view the Saved Graph from the Test parameter screen before re-running the same test.




On selecting GRAPH option, it will display following screen.

In Calibrator Mode, after asking “LOAD NEW BLANK? YES / NO” you will get the following question “LOAD NEW CALIBRATORS? YES / NO”. If you are loading new calibrators and want the results of samples as per new graph, select “LOAD NEW CALIBRATORS? YES”.

Or else, if you are not loading new calibrators and want to use previously stored graph select “LOAD NEW CALIBRATORS? NO”.

If you want to load only ‘Blank’, select “LOAD NEW BLANK? YES” and “LOAD NEW CALIBRATORS? NO”.

In case if you are using previously Stored “Blank” and “Calibrators” then it will print following strings in



**B) In “CUT OFF” Mode:**

NAME TSH		MODE COFF		PRI 450		SEC 0		BLANK 1		NC 3		PC 2		
SMPL 5		RUN		PID		GREYZONE 10.00				CUT-OFF ABSORBANCE 0.160				
	W1	W2	W3	W4	W5	W6	W7	W8	W9	W10	W11	W12		
A	B1 0.050	S3 0.248												
B	NC1 0.006	S4 0.418												
C	NC2 0.012	S5 0.399												
D	NC3 0.115													
E	PC1 0.113													
F	PC2 0.119													
G	S1 0.226													
H	S2 0.219													
											SIMUL		MODIFY	
ESC		SAVE RESULT		LOAD NEXT		PRINT MATRX		PRINT RESULT		ACCEPT TEST				

**When you want to run such an accepted test, select that particular test for instant ‘CUTOFF’ mode and point to ‘RUN’ option present on the run screen.**


This time, the blank and controls will not be loaded directly and you will get the following questions: “LOAD NEW BLANK? YES / NO” and “LOAD NEW CONTROLS? YES / NO”. If you are running both, select ‘YES’, and if you want to use earlier stored data of controls select ‘NO’, so that you can run only samples to get their results without loading Blank and Controls. If you want to load only ‘Blank’, select “LOAD NEW BLANK? YES” and “LOAD NEW CONTROLS? NO”.

In case if you are using previously Stored “Blank” and “Controls” then it will print following strings in the printout-  
 “Using Blank Stored value” and “Using Controls Stored value”

### 11. SAMPLE AND SAMPLE DUPLICATE

While loading the samples you can load a single sample in single well or a single sample in adjacent two well and finally the instrument will take the average of it known as Sample Duplicate. The instrument will load the controls / calibrators automatically in any new test or not accepted test. You have to provide the number of samples. If you are loading single samples keep “DUP”-No. Select “No. of Samples:” to enter total samples.

A. Select touch zone of “SMPL” (Sample)



B. For Duplicate sample


If you want to load duplicate samples, select “DUP N” option before entering number of samples option present in sample screen which shown in below screen.

The screenshot shows a software interface for sample management. At the top, there are several data fields: NAME (TSH), MODE (MSTD), PRI (450), SEC (0), BLK (Y), and CALI (2). Below these are fields for SMPL (0), RUN, and PID, along with a RANGE (10.30-16.61). The main area is a grid with columns labeled W1 through W12 and rows labeled A through H. The first three rows (A, B, C) contain sample identifiers and values: B1 0.000, C1 0.000, and C2 0.000. At the bottom, there are control buttons: ESC, SAVE RESULT, LOAD NEXT, PRINT MATRX, PRINT RESULT, ACCEPT TEST, and GRAPH. A pop-up menu is open, showing 'SMPL 0' and 'DUP N' as the selected options, with 'TABLE', 'SIMUL', and 'MODIFY' also visible.

After selecting “DUP N” it toggles to “DUP Y” and then select “SMPL”.

This screenshot is identical to the one above, but the pop-up menu now shows 'DUP Y' as the selected option instead of 'DUP N'. All other elements, including the data fields, grid, and control buttons, remain the same.

C. Now, enter the no. of sample by Selecting touch zone of “SMPL” (Sample), it will display the numeric screen. After entering no. of sample select “ENT” and then sample will automatically load.



D. After selecting DUP option and entering no. of sample, it displays following screen.

NAME TSH		MODE MSTD		PRI 450		SEC 0		BLK Y		CALI 2	
SMPL 10		RUN		PID		RANGE 10.30-16.61					
	W1	W2	W3	W4	W5	W6	W7	W8	W9	W10	W12
A	B1 0.000	S3 0.000	S7 0.000								
B	C1 0.000	S8 0.000	S8 0.000								
C	C2 0.000	S4 0.000	S8 0.000								
D	S1 0.000	S5 0.000	S9 0.000								
E	S1 0.000	S5 0.000	S10 0.000								
F	S2 0.000	S6 0.000	S10 0.000								
G	S2 0.000	S6 0.000									
H	S3 0.000	S7 0.000									
								TABLE		SIMUL	
ESC		SAVE RESULT		LOAD NEXT		PRINT MATRX		PRINT RESULT		ACCEPT TEST	
								GRAPH			


(\* Similarly, user can load the samples and Sample Duplicates in all other modes)

## 12. PATIENT IDENTITY / SAMPLE IDENTITY ENTRY

To enter Patient ID for all samples, user first run the test. After running test select sample well to assign PID to the sample before “SAVE RESULT”. User can assign the PID to the sample after “ACCEPT TEST” but before “SAVE RESULT”. After entering PID of sample select “SAVE RESULT” option to save the sample result with PID.

NAME TSH		MODE COFF		PRI 450		SEC 0		BLANK 1		NC 3		PC 2	
SMPL 5		RUN		PID		GREYZONE 10.00				CUT-OFF ABSORBANCE 0.160			
	W1	W2	W3	W4	W5	W6	W7	W8	W9	W10	W11	W12	
A	B1 0.050												
B	NC1 0.006	S4 0.418											
C	NC2 0.012	S5 0.399											
D	NC3 0.115	S6 0.000											
E	PC1 0.113	S7 0.000											
F	PC2 0.119	S8 0.000											
G	S1 0.226	S9 0.000											
H	S2 0.219	S10 0.000											
TAG: S3		ESC		TABLE				SIMUL		MODIFY			
ABS: 0.001													
RES: 1.854													
PID: 3													
REM: POS				LOAD NEXT		PRINT MATRX		PRINT RESULT		ACCEPT TEST			

Select PID, it will display alphanumeric screen.



One by one you can enter name of each patient. The character length of PID is only 7 characters.

(\*NOTE: After selecting “SAVE RESULT”, PID get stored in instrument memory. To change the PID of sample you have to select “LOAD NEXT” option)



There is another way to entered the PID which shown in below.

To enter Patient ID for all samples there is an option provided in loading screen shown as “PID”

After entering number of samples either single or duplicate, select “PID”. It will display confirmation message i.e. “Delete Previous PIDs.. YES or NO”, Select YES

A) Individual PID entry -


NAME TSH		MODE MSTD		PRI 450		SEC 0		BLK Y		CALI 2		
SMPL 10		RUN		PID		RANGE 10.30-16.61						
	W1	W2	W3	W4	W5	W6	W7	W8	W9	W10	W11	W12
A	B1 0.000	S3 0.000	S7 0.000									
B	C1 0.000	S8 0.000	S8 0.000									
C	C2 0.000	S4 0.000	S8 0.000									
D	S1 0.000	S5 0.000	S9 0.000									
E	S1 0.000	S5 0.000	S9 0.000									
F	S2 0.000	S6 0.000	S10 0.000									
G	S2 0.000	S6 0.000	S10 0.000									
H	S3 0.000	S7 0.000										
										SIMUL		MODIFY
ESC		SAVE RESULT		LOAD NEXT		PRINT MATRX		PRINT RESULT		ACCEPT TEST		GRAPH


Delete Previous PIDs..

YES

NO

Sample	PID	
S01		ESC
S02		
S03		NEXT
S04		
S05		PRV
S06		
S07		
S08		SEQL





One by one you can enter name of each patient. The character length of PID is only 7 characters in the following combinations –

Since PID of the Patient is entered with the help of Alpha-numeric characters, the combination must be as mentioned below.

Alphabets	Numbers
4	3
3	4
2	5

**For example -**

1. In the first combination, user can enter PID with 4 Alphabets and 3 Numbers. That means, entered PID should be PQRS123
2. Similarly, for the second combination, user can enter PID with 3 Alphabets and 4 Numbers. And therefore, entered PID will be UTV1234
3. And for the third combination, user can enter PID with 2 Alphabets and 5 Numbers ie. PID will be displayed as RB12345
4. But one can't enter the PID with 1 Alphabet and 6 Numbers since its maximum numeric count is 65500.

(Note: The same combinations should be used during SEQUENTIAL entry also.)


### B) Sequential Patient ID Entry

This one by one entry will take more time, therefore one more option to enter Patient Ids is provided. This option is “SEQL” sequential entry.

If you select “PID” option for patient ID entry, the following question appears on the screen:

“Delete Previous PIDs? Yes / No”. If you select NO, it will directly load the previously entered PIDs for sample and if you select YES, it will delete the previous PIDs.

Now select ‘SEQL’, it will display Alphanumeric screen where user can enter different sequential PIDs. Follow the steps below-



Final screen will appear which displays the sequential IDs of patient samples for the well location that is selected.

Sample	PID	
S01	PRST123	ESC
S02	PRST124	
S03	PRST125	NEXT
S04	PRST126	
S05	PRST127	PRV
S06	PRST128	
S07	PRST129	SEQL
S08	PRST130	

Go to the next and previous screen with the help of “NEXT” and “PREV” option present on the screen.

(\*NOTE: The entered PID get stored in instrument memory, until you delete it by using option, “Delete previous PIDs? Yes / NO”)

### 13. SAVE PLATE AND SAVE RESULTS

In all modes of operation, after completing RUN to save the result first select “ACCEPT TEST” and then select “SAVE RESULT” option. “Save result” provide two option i.e. “SAVE RESULT” & “SEND RESULT”.

**A) In Cut Off mode:**

After running the cutoff mode select “ACCEPT TEST” option, it will display the message i.e. “Test Accepted”.

	W1	W2	W3	W4	W5	W6	W7	W8	W9	W10	W11	W12
A	B1 0.050	S3 0.248										
B	NC1 0.006	S4 0.418										
C	NC2 0.012	S5 0.399										
D	NC3 0.115											
E	PC1 0.113											
F	PC2 0.119											
G	S1 0.226											
H	S2 0.219											


Test Accepted...

After test accepted select save result option present in run screen. It will display two option i.e. “SAVE RESULT” & “SEND RESULT”, in that select “SAVE RESULT” option

The screenshot displays the run screen interface with the following data and controls:

NAME TSH		MODE COFF	PRI 450	SEC 0	BLANK 1	NC 3	PC 2				
SMPL 5	RUN	PID	GREYZONE 10.00			CUT-OFF ABSORBANCE 0.160					
W1	W2	W3	W4	W5	W6	W7	W8	W9	W10	W11	W12
A	B1 0.050	S3 0.248									
B	NC1 0.006	S4 0.418									
C	NC2 0.012	S5 0.399									
D	NC3 0.115										
E	PC1 0.113										
F	PC2 0.119										
G	S1 0.226										
H	S2 0.219										

Control buttons at the bottom of the screen include: ESC, SAVE RESULT, SEND RESULT, LOAD NEXT, PRINT MATRX, PRINT RESULT, ACCEPT TEST, TABLE, SIMUL, and MODIFY. A callout box highlights the 'SAVE RESULT' and 'SEND RESULT' options.



S

The screenshot shows the main interface with a 'Result Saved' message in the center. The interface includes a header with parameters: NAME (TSH), MODE (COFF), PRI (450), SEC (0), BLANK (1), NC (3), and PC (2). Below this are fields for SMPL (5), RUN, PID, GREYZONE (10.00), and CUT-OFF ABSORBANCE (0.160). A grid of 12 wells (W1-W12) is shown with rows A-H. Well A1 contains B1 (0.050) and S3 (0.248). Well B1 contains NC1 (0.006) and S4 (0.418). Well C1 contains NC2 (0.012) and S5 (0.399). Well D1 contains NC3 (0.115). Well E1 contains PC1 (0.113). Well F1 contains PC2 (0.119). Well G1 contains S1 (0.226). Well H1 contains S2 (0.219). At the bottom, there are buttons for ESC, SAVE RESULT, LOAD NEXT, PRINT MATRX, PRINT RESULT, and ACCEPT TEST. A 'Result Saved' message is displayed in the center, with buttons for SAVE RESULT, SEND RESULT, TABLE, SIMUL, and MODIFY.

After completing plate and result save process it will display following screen.


This screenshot shows the same interface as the previous one, but with a different set of control buttons at the bottom. The 'Result Saved' message is no longer present. The buttons at the bottom are: ESC, SAVE RESULT, LOAD NEXT, PRINT MATRX, PRINT RESULT, and ACCEPT TEST. The buttons for TABLE, SIMUL, and MODIFY are also present in the lower right area.

1. "SEND RESULT" is used to send the data from instrument to computer through USB or serial RS232 connection.
2. "SAVE RESULTS" option is used to save the details related to samples and it helps the user to save the entire plate information along with the control/calibrator information (maximum 100 plate data can be saved).




## 14. VIEW STORED DATA

The stored plate or stored results appear on the screen by selecting option “View Plate” in “MAIN MENU”.



### 14.1. PLATE RESULT:

14.1.a. It displays the details of stored plate (last 100 plates) with date, test name and mode of operation in descending order. User can select any test to view its plate results by selecting this option.



S

NAME SIN		MODE MSTD													
		W1	W2	W3	W4	W5	W6	W7	W8	W9	W10	W11	W12		
A	C1 0.083														
B	S1 0.544														
C	S2 0.367														
D	S3 0.300														
E	S4 0.700														
F	S5 0.326														
G	S6 0.491														
H															

TABLE

ESC    SEND RESULT    PRINT MATRX    PRINT RESULT    GRAPH


“PRINT MATRX” will print the details of entire plate in matrix format as explained earlier.

“PRINT” will print only the data displayed on screen.

“SEND RESULT” can be used to send the data from instrument to computer through USB or serial RS232 connection.

#### 14.1.b. TABLE:

It shows the result in tabular form. It displays the details of stored plate (last 100 plates) with well no., Tag, unit, PID, remark, absorbance and there result.



S

NAME TSH		MODE COFF		PRI 450		SEC 0		BLANK 1		NC 3		PC 2	
SMPL 5		RUN		PID		GREYZONE 10.00				CUT-OFF ABSORBANCE 0.160			
	W1	W2	W3	W4	W5	W6	W7	W8	W9	W10	W11	W12	
A	B1 0.050	S3 0.248											
B	NC1 0.006	S4 0.418											
C	NC2 0.012	S5 0.399											
D	NC3 0.115												
E	PC1 0.113												
F	PC2 0.119												
G	S1 0.226												
H	S2 0.219												

TABLE      SIMUL      MODIFY

ESC      SAVE RESULT      LOAD NEXT      PRINT MATRX      PRINT RESULT      ACCEPT TEST

NAME TSH		MODE COFF				
WELL_NO	TAG	ABS	RESULT	UNIT	PID	REMARK
A01	B1	0.050	10.4		VSR1	EQU
B01	NC1	0.006	12.6		VSR2	HI
C01	NC2	0.012	7.450		VSR3	LO
DO1	NC3	0.115	13.80		VSR4	HI
E01	PC1	0.113	267.7		NRA1	EQU
F01	PC2	0.119	220.3		NRA2	LO
G01	S1	0.226	180.2	mg/dl	NRA3	LO
H01	S2	0.219	326.7	mg/dl	NRA4	HI


PRV      NXT      MATRIX

ESC      SEND RESULT      PRINT MATRX      PRINT RESULT

## 14.2. PATIENT RESULTS

It will show the details of all samples, saved in memory using option “SAVE RESULTS” in all different modes, as explain earlier. This option will display the data in table format as shown below. Test name, Mode of test operation, result of sample, PID and Remark will appear on the screen. The instrument has a memory to store 2500 sample results.

This option doesn't show the information of controls or calibrators.



“PRN” is used to print the details as available on the screen.

“NXT” and “PRV” options are used to browse through the pages so that one can view all the results one after the other.

## NOTE:

1. It's display alert message for last 10 plate data to store the results.  
i.e. "Only 10 plates remaining"
2. After completing 100 plate data it will display the message i.e.  
"Record full"
3. When there is not enough memory the following message appears:  
"First plate data Deleted"

The instrument automatically deletes a number of samples that are stored in the beginning and then save the new samples at that location.


4. If a test is being deleted or edited then previously stored results of that particular test are displayed with their test mode as "X". But test name will display as 'Deleted'.

16 / 08 / 16					
TEST_NAME	MODE	RESULT	UNIT	PID	REMARK
DELETED	X	10.48		VSR1	EQU
DELETED	X	12.63		VSR2	HI
DELETED	X	7.450		VSR3	LO
DELETED	X	13.80		VSR4	HI
SIN	S	267.7		NRA1	EQU
SIN	S	220.3		NRA2	LO
SIN	S	180.2		NRA3	LO
SIN	S	326.7		NRA4	HI

ESC    NXT    PRV    PRN

## 15. UTILITIES


UTILITIES consists of 10 different functions which can be viewed as follows-



### 15.1. Date & Time Setting

This Utility helps the user to set current Date and Time.

Set current Date and Time by Selecting DD, MM, YY and Hrs, Min, Sec.



To enter /change the date and time, touch on the shaded zones accordingly.


A numeric submenu screen will be displayed.

Enter the current date and time.

On completion, Touch ESC to escape


### 15.2. Program Filter

It consists of two editable filters and thus helps the user to edit the filters if any.



### 15.3. Printer On / Off setting:


This option is use for multi-purpose. Using this option user can use internal printer, external printer or user can switch-off the printer






### 15.4. Enter Clinic Name:


Enter your Clinic Name and select “ENTER”. The character length is maximum 20 characters.



### 15.5. Shaker (Plate Shaking Mode)

Shaker is used for shaking the plate before running any test. The Shaker screen is displayed as follows –






The total shake time can be entered in seconds

Regarding the “Shake Speed”, there are 10 types of speeds. You can choose any speed by entering numbers from 1 to 10. The shaking speed increases as the number increases.


Speed	Respective RPM
1	200
2	210
3	220
4	230
5	240
6	250
7	260
8	270
9	280
10	290

After entering shake time and shake speed, select Start option. When “Start” is selected, the tray will come out and on the screen you will obtain a message “Is plate loaded? Yes / No”. Load the plate and select YES option.



After selecting ‘YES’ the tray will go inside and instrument will do the shaking of plate for the given time period with selected speed.


When the shaking is finished, the plate holder will come out and display a message- “Is plate removed? Yes / No”. Remove the plate and select ‘YES’ to end the process.



Similarly, to stop the Shaking process in between select “STOP” option present on the screen. After selecting STOP option it will display the message string – ‘Is Plate Removed? YES / NO’. Select NO to terminate the ongoing process.

### 15.6. Communication:

It is a toggle between USB and SERIAL communication. User can either switch to USB or Serial communication by connecting their respective cables.



This setting is very important, whenever you want to transfer data from instrument to computer. You can transfer data either using USB or by Serial RS232 cable.

### 15.7. CV Diagnostics: OFF

**NOTE:** The option “CV DIAGNOSTIC OFF” is not for users. It is for factory use only.

### 15.8. UPDATE FIRMWARE


**NOTE:** The option “UPDATE FIRMWARE” is not for users. It is for factory use only.

### 15.9. BI-DIRECT:

This software use to transfer the data for LIS software.


**15.10. PASSWORD:**

Select Password option in order to Set the Password so that one can't ADD, EDIT, DELETE and SEND / RECEIVE Test records without entering the valid Password.




**A) ENABLE / DISABLE:**

Enable / Disable helps the user to set the Password or Remove the Password from the Password Utility. Enter Password with the help Alphanumeric screen which is displayed after selecting Enable / Disable option.




**B) SET PASSWORD:**

This option helps to set the PASSWORD along with the USER name.




**C) LOGIN:**

After entering the Password user needs to select LOGIN.



To Login/ enable user, enter the USER name and select ENTER


**D) LOGOUT:**



To Logout / disable user, select LOGOUT option present on the screen

## 16. COMMUNICATION

The instrument is equipped with an RS232 serial port for PC configuration ie. User-computer interface. A cable is available to link the instrument to PC.




Communication will only start when both ends detect the presence of an active terminal or device.

RS232 port settings in a windows Operating system.

PORT SETTINGS	
Bits per second	: 9600
Data Bits	: 8
Parity None	: None
Stop Bits	: 1
Flow control	: None

There is an option present on the Main menu screen named as COMMUNICATION which is used for transferring data from instrument to the computer and also from instrument to instrument.






**16.1. Send Results:**

**SEND RESULTS** helps the user to transfer all the Results stored in the instrument to the computer via Serial Communication.

Results can be send either through USB or SERIAL communication.

Before Sending the Results check whether the Serial cable / USB cable is connected to both the ends of the instrument and the computer respectively.




These Results are send through Application which are created in “.csv” format, Convert the corresponding format to Excel file which is created as shown below –

The screenshot shows a Microsoft Excel spreadsheet titled "Microsoft Excel - PRL0G2012831\_1". The spreadsheet contains a table of test results. The columns are labeled A through H, with the following headers: A: Date, B: TestName, C: Mode, D: Result, E: Unit, F: PID, G: Remark:-. The data rows show test results for various dates (all 31/08/12) and test names (all SIN). The results include values for Mode (S), Result (ranging from 18.93 to 113.5), Unit, and PID (CID1 to CID14). The Remark:- column is empty for all entries.


	A	B	C	D	E	F	G	H
1								
2	TEST RESULTS:							
3								
4	Date	TestName	Mode	Result	Unit	PID	Remark:-	
5								
6	31/08/12	SIN	S	113.5		CID1		
7								
8	31/08/12	SIN	S	116.6		CID2		
9								
10	31/08/12	SIN	S	96.42		CID3		
11								
12	31/08/12	SIN	S	78.16		CID4		
13								
14	31/08/12	SIN	S	18.93		CID5		
15								
16	31/08/12	SIN	S	78.84		CID6		
17								
18	31/08/12	SIN	S	67.02		CID7		
19								
20	31/08/12	SIN	S	111.9		CID8		
21								
22	31/08/12	SIN	S	101.9		CID9		
23								
24	31/08/12	SIN	S	90.22		CID10		
25								
26	31/08/12	SIN	S	64.82		CID11		
27								
28	31/08/12	SIN	S	34.25		CID12		
29								
30	31/08/12	SIN	S	73.61		CID13		
31								
32	31/08/12	SIN	S	72.27		CID14		
33								

### 16.2. Send Test Records:

This option helps the user to send the Test Records created in one instrument to another instrument serially with the help of serial cable. And therefore one has to select serial communication in Utilities.




Suppose by default if Communication is set to USB then it will display following screen with msg string.



### 16.3. Receive Test Records:

This option is used for receiving the Test Records from one instrument to another. On selecting “RECV TEST RECORDS”; it will display the string – Receive Test Records? Yes / No. Select “Yes” in order to receive the data from the other instrument.




Note: 1. User must first restart the instrument before sending the data from one instrument to another.

2. Remember that while transferring Test Records from one instrument to another one must first select “RECV TEST RECORDS” and then select “SEND TEST RECORDS” in order to send the data from one instrument to another.

For instance, let us take an example which will show how the Test Records are transferred from instrument to instrument.


**16.4. Check Serial Communication:**



**This option is only applicable for Serial Communication. If Communication is set to USB then it will display the message i.e. "Serial comm problem!!!!"**

## 17. TECHNICAL DIAGNOSIS

Technical Diagnosis is present in “Main Menu”.




### 17.1. Clear Memory

This option is used to clear entire saved test records and entire saved patient results along with saved plates. The “Test Records” are the different kind of programmed test that you have made in different modes of operation and saved in instrument program memory.

When you select Clear Patient Results option, the instrument will first ask, “Clear Results? YES / NO”. If you select “YES” it will start clearing only saved patient results and saved plates, stored in instrument memory, and shows message “Please wait..”.

This option does not clear the programmed test saved in instrument memory.


When you select Clear Test Memory option, this time instrument will ask, “Clear Tests? YES / NO”. If you select “YES” it will start clearing all programmed test saved in instrument memory, and shows message “Please wait ..”.



### 17.2. Lamp Amplitude


This option is used to check the light gain of filters used in the instrument. There are 2- channels used of same type. This option shows the gain of each channel filter in table format. In the instrument the range of amplitude is from 800.00 to 2450.00. If the amplitude goes below 800.00 for any channel, the message “Amp –Low” will appear on the screen for that particular channel if the amplitude goes more than 2490.00, than the over range saturation voltage figure of 2500.00 will be passed.

The gains of filters are set within a required range. When you select this option, the lamps will turn ON, tray holder plate will comes out and the instrument will start showing amplitudes (gains) of all filters one by one.





2. Filter Amplitude:



**17.3. Printer Self test**

It is used to check the printing of thermal printer. Is it printing perfectly or not?  
When you select this, it will print the message i.e. 'printer check ok'.

**17.4. Lamp Centering Check:**

**NOTE: "Lamp Centering Check" is not for users. It is only for service purpose.**

**17.5. PLATE HOME:**

**NOTE: "PLATE HOME" is not for users. It is only for service purpose.**

**17.6. PLATE PARK:**

**NOTE: "PLATE PARK" is not for users. It is only for service purpose.**

**17.7. FILTER HOME:**

**NOTE: "FILTER HOME" is not for users. It is only for service purpose.**

**17.8. LAMP HOME:**

**NOTE: "LAMP HOME" is not for users. It is only for service**

**17.9. CHECK MOTOR MOVEMENT:**

**NOTE: "CHECK MOTOR MOVEMENT" is not for users. It is only for service**

**17.10. FILTER GRAPH:**

**NOTE: "FILTER GRAPH" is not for users. It is only for service**

## 18. Trouble Shooting

	CAUSE /CORRECTIVE ACTION
<p>1) Printer disabled Disable Printer YES / NO?</p>	<p>You will get this message, if paper is not loaded properly or lever is not at correct position. Check all possibilities. Also check ON LINE and FEED LED glowing or not and when you switch 'ON' the instrument check paper is moving forward or not. If LEDs are not glowing and paper is not moving forward contact service engineer.</p>
<p>2) "Check Light Path!!! CH x"  Where 'x' is any channel number X = 1/2</p>	<p>You can get this message in any mode of operation. Before reading absorbance, instrument is checking the filter reference voltages of all 2 channels. If filter voltage of any channel is less than minimum required voltage, it will display a message "Check Light Path!!! CH x". Here 'x' is channel numbers having reference voltage less than minimum required.</p> <p>For ex. If reference voltage of channel 2 is less the message will be displayed "Check Light Path!!! CH 2". It will indicate all channel numbers having less reference voltage, like "Check Light Path!!! CH 1 2 "</p> <p>In such a case, there is a possibility of filter gain of that channel is reduced or intensity of lamp of that channel has become poor, so contact factory engineer.</p>
<p>3) Invalid Assay</p>	<p>If in the Cut Off and Multi standard mode the controls and standards are not ok then this error message will appear.</p>
<p>4) Memory Full</p>	<p>If Number of saved tests exceeds the memory limit of 250 tests then delete the unwanted tests and save the test.</p>
<p>5) "Check USB Application on Computer</p>	<p>Using data receiving application on computer, 1) If you have not connected the data cable (USB or Serial) correctly, you will get this message. 2) If you have selected 'Communication Setting', USB on instrument, and you have selected a 'Serial' communication on computer data receiving application, you will get this message. In such a case select 'USB'</p>

<p>Do you want to continue? YES / NO”</p>	<p>communication, on computer and select ‘YES’ key to transfer data. To avoid this message, connect proper data cable, select the same communication medium (USB or Serial), on computer data receiving application, which you have selected on instrument ‘Communication Setting’.</p>
<p>6) Filter Movement Error!!! OR Filter Index Error!!!</p>	<p>Instrument will display this message, if there is problem with MOC sensor PCB or with filter tray stepper motor. The MOC sensors are used to detect the position of filter tray. These two sensors are mounted on bottom of the mechanism. One detects the home position and second detects the each filter position of filter tray. In such a condition we need to check the connection to these sensor and also the connection of stepper motor, which is used to move filter tray.</p>
<p>7) Plate Movement Error!!! OR Plate Index Error!!!</p>	<p>Instrument will display this message, if there is problem with MOC sensor PCB or with plate carrier stepper motor. The MOC sensors are used to detect the position of plate carrier. These two sensors are mounted on top of the mechanism. One detects the home position and second detects the each strip position of plate carrier. In such a condition we need to check the connection to these sensor and also the connection of stepper motor, which is used to move plate carrier.</p>
<p>8) Please Enter Filter value</p>	<p>Instrument will display this message, if primary filter not entered</p>
<p>9) Error-Same Filter</p>	<p>Instrument will display this message, if there is Same filter</p>
<p>10) CheckSum Error</p>	<p>Instrument will display this message, if sum is not proper user need to resend the command.</p>

---

## 19. DECONTAMINATION

### 19.1. Decontamination Procedure

---

If the instrument is to be shipped after being exposed to potentially hazardous material, it should be decontaminated. The following procedure outlines the method of decontaminating the instrument before packaging and shipment.

### 19.2. Purpose of Decontamination

---

Decontamination minimizes the risk to all who come in contact with the instrument during shipping, handling, and servicing.

### 19.3. General Considerations

---

- Any laboratory instrument that has been used for clinical analysis is considered a biohazard and should be decontaminated prior to handling. Intact skin is generally considered an effective barrier against infectious Organisms; however, small abrasions and cuts may not be always being visible. Prophylactic gloves must be worn when handling instruments that have not been decontaminated. Gloved hands should be considered contaminated at all times and must be kept away from eyes, mouth and nose at all times.
- Mucous membranes are considered prime entry routes for infectious agents. Wear eye protection and a surgical mask when there is a possibility of aerosols.
- Eating and drinking while decontaminating instruments is not advisable.

### 19.4. Procedure

---

- A solution of .5% Sodium Hypo Chlorite (NaOCL) solution (Bleach) is used. Commercial bleach is 5% NaOCL; household bleach is 3% NaOCL. When using commercial bleach, use a 10:1 mixture; if using household bleach, a 6:1 mixture is required. This is a caustic solution. It is important to wear gloves and eye protection when handling it.
- Wipe down the carrier and all exposed surfaces of the unit with the bleach solution. Remove the top shroud of the instrument and wipe down the top surface of the instrument base, as well as the inside of the top shroud.
- Reassemble the unit and discard the used gloves and towels.

---

## 20. SAFETY CLEARANCE CERTIFICATE

Please complete all information requests on this form prior to returning the instrument to the manufacturer or your local distributor for servicing, repairs or return. Thank you for your co-operation.

Customer \_\_\_\_\_ Contact \_\_\_\_\_

Address \_\_\_\_\_ Position \_\_\_\_\_

Dept \_\_\_\_\_

Tel: \_\_\_\_\_

Country \_\_\_\_\_ Fax: \_\_\_\_\_

Post Code \_\_\_\_\_

Model No. \_\_\_\_\_ Serial No. \_\_\_\_\_

Accessories \_\_\_\_\_

Returned \_\_\_\_\_

Date of Purchase (if known) \_\_\_\_\_

Complaint \_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

---

Has the equipment been exposed to any of the following: (\*delete as applicable)

a) Blood, body fluids, pathological specimens \*YES/NO

If YES, please specify \_\_\_\_\_

\_\_\_\_\_

b) Other Biohazard \*YES/NO

If YES, Please specify

\_\_\_\_\_

\_\_\_\_\_