



MICRO SHAKE 1200

ELISA Plate Shaker



Synchronised Harmony with Skaking & Incubation

FEATURES :

- Simultaneous orbital shaking & incubation operation
- Buzzer indication on completion of incubation
- Indication of Remaining Time
- Display of temperature on pressing TEMP Key



RPM
400 to 700



Incubator



Shaker

TECHNICAL SPECIFICATIONS :

Operating Modes	:	Shaker Incubator Shaker & Incubator
Temperature Control		
Temperature Range	:	370 C to 420 C
Resolution	:	10 C
Incubation Time	:	1 to 999 min
Shaker		
Frequency	:	400 to 700 RPM
Operating Position	:	On horizontal flat, rigid and vibration free surface
Operating Conditions		
Temperature	:	From + 180 C to 320 C
Relative Humidity	:	Up to 80%
Storage Conditions		
Temperature	:	From - 100 C to 350 C
Relative Humidity	:	Up to 80%
Enclosure	:	ABS Fire retardant
Size (cm)	:	28 X 25 X 16 (l X b X h)
Weight	:	3 Kg (Approx)

Due to continuous development, specifications are subject to change without prior notice.

Global Diagnostics B

Sijsjesstraat 4 2440 Geel, BELGIUM, Tel.: +32468220039, info@globaldiagnosticsb.com